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

- EN**
1. Centrifugal Blower
  2. Exhaust flexible tube
  3. Laminar flow chamber or plenum
  4. Laminar flow HEPA filter
  5. Exhaust HEPA filter
  6. Working area
  7. Outside air inlet and protection barrier
  8. Exhaust outlet

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
**EN Annex A3. CONTROL ALARM**

1. On/Off button
2. Navigation buttons
3. Graphic screen

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


**EN**

- |                    |                               |
|--------------------|-------------------------------|
| 1. Connecting pipe | 3. Active charcoal filter box |
| 2. Support         | 4. Safety Cabinet             |

- A. Place the connecting pipe (1) on top of the safety cabinet air extraction. It will be sealed by the adhesive tape provided.
- B. Place support (2) on the wall nearest the safety cabinet and hang up the active charcoal filter box (3). This can also be placed in any available position.
- C. Apply the required amount of suitable diameter flexible piping.

**IMPORTANT:** It is mandatory to install an additional blower to guarantee a proper exhaust flow

**GS has not been tested for Active Carbon Filter**

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**EN**

*For worktable and sides, use sterile or non-use disposable cloths, which do not shed particles or fibres.*

*For disinfection, these cloths will be used, slightly dampened with a non-corrosive disinfectant solution (70% alcohol solutions), that do not damage or affect the stainless steel and the glass. Cleaning can also be carried out with water and soap prior to disinfection.*

*To disinfect correctly, the appropriate disinfectant solution should be used for the type of work that is being done. Many of these solutions require certain application times depending on the type of micro-biological agent used in the cabinet; we recommend that you consult the instructions for application and precautions for use given by the manufacturer of the disinfectant.*


**Fig. 1**

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**EN**

*To accede to corners or difficult access areas, another type of tools will be used as shown in the Fig 2.*

**Fig. 2**

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**EN**

*The perforated stainless steel worktable is divided into three parts to facilitate cleaning of the interior container and provide for sterilisation of these parts into a vessel. Remember that every segment or shelve has to be cleaned and/or disinfected both in its lower and upper sides. The parts will be cleaned inside the cabinet while this is running.*

**Fig. 3**

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**EN**

*Also be sure of carrying out the cleaning/disinfection in the lower part of the working area (surface and corners).*

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**EN**

*Inside, safety cabinet cleaning is recommended on the following occasions:*

- a) *Before starting any work in the safety cabinet.*
  - b) *After working in the safety cabinet.*
  - c) *Whenever there is a change of work programme.*
  - d) *In event of liquid spilling on the work table.*
  - e) *Before carrying out mechanical or biological control tests in the working area.*
-

**FR Annexe A5.1. NETTOYAGE ET  
DÉSINFECTION SUPERFICIELLE****DE Anhang A5.1. REINIGUNG UND  
DESINFEKTION DER OBERFLÄCHEN****FR**

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**EN**

*For the cleaning of front glass you must use sterile textures or disposable ones without particle or fibres loosing.*

*For disinfection, use the previous textiles slightly dampened with a non-corrosive disinfectant solution (70% alcohol solutions), that do not damage or alter the glass. Cleaning can also be carried out with water and soap prior to disinfection*

*To disinfect correctly, the appropriate disinfectant solution should be used for the type of work that is being done. Many of these solutions require certain application times depending on the type of micro-biological agent used in the cabin; we recommend that you consult the instructions for application and precautions for use given by the manufacturer of the disinfectant.*

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**EN**

*For the glass exterior part cleaning, it is recommended to leave the glass in closed working position. In this position the major part of the glass is accessible. For the interior part of the glass cleaning, it is recommended to leave the glass in its working position in order to assure the user's protection.*

**WARNING** 

*It's important to clean and disinfect the cabinet's working area, with the glass in working position and the cabinet in performance. If it's necessary to separate the frontal-glass mechanism to clean the areas of difficult access, bring the glass to its closed position to set free the frontal mechanism. By separating the frontal-glass, the inner surface of the glass is totally accessible. There must be taken the necessary precautions because the user's protection is lower with the glass separated from cabinet.*

*It's necessary that the glass is kept in closed position before moving back the frontal-glass into the closed position. If not, the glass will not adjust with the frame and user's protection would be seriously affected.*

**WARNING** 

Do not put the hands in the laterals of the frontal support, in order to avoid any entrapment risk.

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
EN

**WARNING** 

**To separate the frontal-glass, pull always only the handles "B" of the metal support, never using the handles of the glass, that are specifically designed for the vertical movement of the glass.**

**WARNING** 

The glass is required to remain in the closed position before closing the frontal towards the cabinet. In case this is not fulfilled, the glass will not be properly adjusted with the frame, and the user protection will be seriously affected.

<b>FR</b>	<b>Annexe A5.2. NETTOYAGE</b>	<b>ET</b>	<b>DE</b>	<b>Anhang A5.2. REINIGUNG</b>	<b>UND</b>
	<b>DÉSINFECTION VITRE FRONTALE</b>			<b>DESINFEKTION DER FRONTSCHIEBE</b>	

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**EN**

*Gas decontamination of a Biological/Cytostatic Safety Cabinet is necessary in the following cases:*

- *Before undertaking any safety cabinet maintenance work.*
- *Before removal of HEPA filters for replacement.*
- *Before carrying out mechanical or biological control tests in the safety cabinet.*
- *Before moving the safety cabinet from one work area to another.*
- *Before changing a work program.*
- *Whenever liquid is spilled on the work table.*

*The usual gas for total decontamination is formaldehyde, made by depolymerisation of paraformaldehyde (solid) or formol (liquid).*

*The quantity of formaldehyde to be used is calculated according to Safety Cabinet total volume (see table 1).*

	<b>Bio II Advance 3</b>	<b>Bio II Advance 4</b>	<b>Bio II Advance 6</b>
<b>Volumen total (m<sup>3</sup>)</b> <b>Total volume (m<sup>3</sup>)</b>	0,75	1,00	1,50
<b>Formaldehído (mg)</b> <b>Formaldehyde (mg)</b>	45	60	90

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**EN**

*Sufficient formaldehyde is required to produce a concentration in the air of at least 50 mg/m<sup>3</sup>.*

*The formaldehyde vapour can be produced by :*

- *Evaporation of formalin (percentage by volume of formalin in an aqueous solution: ± 36%).*
- *Thermal depolymerisation of paraformaldehyde, provided there is an adequate level of humidity.*
- *NOTE: You must spray 60 ml of formalin with 60 ml of water per cubic meter of cabinet volume.*

**WARNING**

*Formaldehyde steam is explosive when the percentage in volume of the mixture with dry air reaches 7.75%; its flash point is 430°C. The national limits on exposure in the workplace should be taken into account.*

*The necessary material to perform decontamination by gas is as follows:*

- *Formaldehyde.*
- *A formolizer/neutralizer unit capable to perform automatically both formolization and neutralization in sequence (such as Telstar unit Formonia VAP-3).*
- *An hermetic plastic cover (PVC or equivalent)*
- *PVC sheet or PVC laminated aluminum sheet*
- *Duct tape*

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**EN****WARNING**

*Formaldehyde penetrates with difficulty and his efficiency depends on the temperature and humidity. It is more effective over 20°C and with relative humidity of 65%. The use of excessive quantities can cause polymer deposits inside the cabinet and can contribute to the filter blockade.*

*It is recommended that workers use proper protection to make decontamination operation, in order to protect them from an accidental formaldehyde gas exposure. Mainly, they must cover the eyes with leak tight safety goggles and, when necessary, use respiratory protection with safety hood. Formaldehyde steams are very irritant for the eyes, throat and nose. Nevertheless, performing decontamination work with enough precautions means that the operator will not face any contact with the gas.*

*Steps to follow:*

- a) *Refer to local regulation concerning the use of formaldehyde.*
- b) *Fumigation must be performed with appropriate precautions and by a qualified person with good knowledge of formaldehyde usage. At any case, a visible warning should indicate the cabinet is being fumigated.*
- c) *Measure the needed quantity of formaldehyde to be used and put it in the formolizer (see table 1).*

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**EN**

- d) *The safety cabinet should be sealed before fumigation to ensure that no formaldehyde can leak into the laboratory or the other rooms. For that, cover the safety cabinet with the plastic cover. This plastic cover is completely airtight and has been tested against leaks with the ammonia test. It consists of testing the plastic cover at a pressure of 80 Pa with an ammonia atmosphere inside. Externally, an ammonia detector checks no leakage. Check that the PVC sheet or PVC laminated aluminum sheet is continuous (one single piece without any joint)*
- e) *Place the PVC sheet or PVC laminated aluminum sheet below the table support of the cabinet (by raising first one side and then the other side)*
- f) *Switch on the cabinet, disable the socket inlets and switch off the fans*
- g) *Separate the frontal of the cabinet and raise the glass to the maximum position*
- h) *Plug in the vaporized with ammonia into the working surface while keeping the plug switched off by the display control*
- i) *Cover totally the cabinet using the hermetic PVC bag and letting go out only the cabinet power cord by the lower part of the bag*

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
**EN**

- j) Seal the cabling together with the PVC bag and the PVC sheet, with duct tape  
 Sticking process details:
- Duct tape with Omega shape covering the cabling 360° (fig. 1 & 2).
  - It is necessary to turn several times the cabling to guarantee a good sealing all around the 360° (fig. 3)
  - Raise the cabling by the free exterior part (fig. 4) and put another piece of duct tape (in red in fig. 5)
  - Stick this piece of duct tape against the former piece of duct tape with the Omega shape (fig. 6) and turn several times the cabling to guarantee a good sealing (fig. 7)

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**EN**

- k) *Complete seal the plastic cover, against the PVC sheet base, with duct tape, taking special care about avoiding any folding that may lead to some possible leak (fig. 10-11-12-13)*
- l) *Cover with a complete second round of duct tape the joint between the plastic bag and the PVC sheet base (fig. 14 & 15)*

Fig 10

Fig 11

Fig 12

Fig 13

Fig 14

Fig 15

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**EN**

- m) Go to the formolizer menú and switch on the sockets by means of the display control*
- n) Go on the formolization procedure. Switch manually the fans after 10 minutes for 15 seconds, once the ammonia is partially vaporized*
- o) Once the ammonia is fully vaporized (about 15 minutes), switch again the fans for 15 seconds and look for leaks using the smell, and in case of not finding using the smell, try with a wet PH paper, and look for a colour alteration*

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**EN**

- p) *For a suitable fumigation, it should be ensured that the downstream side of the HEPA filters and the interior of the cabinet are exposed to formaldehyde enough time to ensure inactivation of microorganisms that have penetrated into the filter. This requires an exposure time of over 6 hours. Fumigation is most conveniently done overnight.*
- q) *Even after the treatment, the HEPA filters can only be considered to be 'safe to handle when using the appropriate protective clothing' instead of no sterile. They should be treated in an autoclave or incinerated on removal. Following the fumigation, the seals of the cabinet should be removed and the extractor fan switched on for a suitable period of time to eliminate the residual formaldehyde steam before using the cabinet again or carrying out any maintenance task.*
- r) *For cabinets that recycle air into the room, the fumigation system used should ensure that no member of personnel is exposed to higher levels than those specified by the relevant authorities. For example, by using a suitable neutralisation system or temporary ducts to extract the exhaust fumes into the exterior. Suitable protective respiratory equipment designed for being used with formaldehyde should be available during fumigation in case of an emergency.*
- s) *Once the process has finished, switch off the equipment in the same way, through the control panel.*

**FR** Annexe A5.3. DÉSINFECTION AU **DE** Anhang A5.3. DESINFEKTION MIT **FORMOL** **FORMALDEHYD**

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

**FR** Annexe A7. SCHEMA  
ELECTRIQUE 230V 50/60Hz


**DE** AnnexA7. SCHALTPLAN  
230V 50/60 Hz


EQUIPMENT: <b>BIO-II-ADVANCE</b>		Date: 13/05/2011 Report of drawings: 1005418 Pag 1	
PROJECT REFERENCE: 5049		BIO-II-ADVANCE	
CUSTOMER: TELSTAR TECHNOLOGIES, S.L.		TELSTAR TECHNOLOGIES, S.L.	
DRAWING TYPE: ELECTRIC			
DRAWINGS: 10			
MATERIAL O.F.:			
PROJECT LEADER: Victor Lazaro			
AUTOMATION LEADER: X. Cañellas			
DRAWN BY: X. Cañellas	DATE: 13/05/2011		
APPROVED: X. Cañellas	DATE: 13/05/2011		
REVIEWED: X. Cañellas	DATE: 13/05/2011		
AS BUILT:	DATE:		
REVISIONS:			
ID: A	DATE: 10/04/12	DRAWN BY: X.Cañellas	
ID:	DATE:	DRAWN BY:	
ID:	DATE:	DRAWN BY:	
ID:	DATE:	DRAWN BY:	
ID:	DATE:	DRAWN BY:	
		PROPERTY OF TELSTAR REPRODUCTION FORBIDDEN	PROJECT NAME: 1005418


Drawing	Description	Creating	Modification
01	DRAWING REPORT	21/01/2010	21/01/2010
02	POWER	21/01/2010	21/01/2010
03	CONTROL BOARD	21/01/2010	21/01/2010
04	MOTORS AND SOCKETS	21/01/2010	21/01/2010
05	LAMPS AND UV	21/01/2010	21/01/2010
06	ELECTROVALVES	21/01/2010	21/01/2010
07	MICROSWITCHES	21/01/2010	13/05/2011
08	PROBE	21/01/2010	10/04/2012
09	DISPLAY	21/01/2010	21/01/2010
10	BILL OF MATERIALS	21/01/2010	10/04/2012




Revision				Date				Draw by				Date				Project description				Station				Zone				Drawings: 6 / 10				Previous: 05				Next: 07			
								13/05/2011				X. Cañellas				BIOII-ADVANCE				TELSTAR TECHNOLOGIES, S.L.				ELECTROVALVES				1005418				5049				06			
								X. Cañellas																															
A				B				C				D				E				F				G				H											

Revision	Date	Draw by	Date	13/05/2011	Project description	BIO-II-ADVANCE	Supplier	Zone	Drawings: 9 / 10	Previous: 08	Next: 10
			Draw by:	X. Carleñas			Customer:	TELSTAR TECHNOLOGIES, S.L.		Drawing number:	1005418
			Approved:	X. Carleñas			Drawing description:	DISPLAY		Project reference:	5049
										Drawing number:	09

Date		13/05/2011		Bill of Materials Page		1	
BIO-II-ADVANCE							
TELSTAR TECHNOLOGIES, S.L.							
Symbol	Drawing	Column	Description	N. Refer.	Catalogue	Reference	
-F1	02	D	Filter PREFILTER	FFH-6Z	ESPE	25183	
-BE1	02	D	Standard chassis plug	RS-481-823	ESPE	11290	
-TC1	03	C	Control board CORTEX-M3		ESPE	65401	
-PS1	03	G	Power Supply MEANWELL	S-60-24	ESPE	65436	
-FAN1	04	C	Fan with capacitor FANS-TECH	FS146F	ESPE	67264	
-FAN2	04	D	Fan with capacitor FANS-TECH	FS146F	ESPE	67264	
-BE2	04	E	Socket plug	B012-H	ESPE	510032	
-BE3	04	F	Socket plug	B012-H	ESPE	510032	
-LT1	05	D	UV KIT		ESPE		
			UV lamp BIO-II-ADV 3				
			UV lamp BIO-II-ADV 4				
			UV lamp BIO-II-ADV 6				
			UV reactance BIO-II-ADV 3				
			UV reactance BIO-II-ADV 4				
			UV reactance BIO-II-ADV 6				
-L2 L3	05	D	LIGHTING		ESPE		
			Lamp BIO-II-ADV 3			36408	
			Lamp BIO-II-ADV 4			71424	
			Lamp BIO-II-ADV 6			72352	
			Ballast BIO-II-ADV 3			71421	
			Ballast BIO-II-ADV 4			71421	
			Ballast BIO-II-ADV 6			72354	
-FC2	05	C	Microswitch OMRON	V-155-1A5	ESPE	71692	
-FC1	05	F	Microswitch OMRON	V-155-1A5	ESPE	71692	
-FC4	02	F	Microswitch OMRON	D2F-L2-D	ESPE	71429	
-FC3	03	F	Microswitch OMRON	D2F-L2-D	ESPE	71429	
-SP1	08	C	Airflow probe ELEKTRONIC	EE576	ESPE	800012	
-DIS1	09	C	Display MERLIN		ESPE	66189	
			Total elements	21			

Revision	Date	Draw by	Date	13/05/2011	Project description	Supplier	Zone	Drawing: 10/10	Previous: 09	Next:
A	10/04/12	X. Cañellas	Draw by	X. Cañellas	BIO-II-ADVANCE	Customer:	TELSTAR TECHNOLOGIES, S.L.	Drawing description:	BILL OF MATERIALS	Drawing number:
		Approved:	X. Cañellas							1005418
								Project reference:	5049	Drawing number:
									10	

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**EN 1. PURPOSE**

Check the integrity of the downflow filters.

**2. SCOPE**

Bio II Advance

**3. TOOLS**

- Aerosol generator (cool or hot).
- ATI PAO-4, Shell Ondina 917 fluid or Emery 3004 oils.
- Photometer (ATI 2H).

**4. METHODOLOGY**

Introduce fume by the generator in the working area (Figure 1)

- Turn on the photometer and adjust it for measurements in accordance with the manufacturer's instructions (Figure 2)

**Fig. 1**

**Fig. 2**

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**EN**

- With the nozzle of the probe approximately 2.5 cm from the surface, scan the downstream side of the HEPA filters, including the perimeter of each filter pack, by passing the photometer probe in slightly overlapping strokes over the entire surface.

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**EN**

- Also scan the entire periphery of the filter and the junction between filter and filter mounting frame at a scanning rate that does not exceed 5 cm/s.

**Fig. 3**

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**EN 5. ACCEPTANCE CRITERIA**

- The photometer signal after the filter should be smaller or equal to a 0,01% of the specific concentration before the filter.

**Figura 4.** Resultado correcto: no fuga  
Correct result: filter not leaking

**Figure 5.** Wrong result: filter leaking

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**EN 1. PURPOSE**

*Prove the not existence of leaks in the HEPA exhaust filters*

**2. SCOPE**

*Bio II Advance*

**3. TOOLS**

- *Aerosol generator (cool or hot).*
- *ATI PAO-4, Shell Ondina 917 fluid or Emery 3004 oils.*
- *Photometer (ATI 2H).*

**4. METHODOLOGY**

*Introduce fume by the generator in the working area (Figure 1)*

- *Turn on the photometer and adjust it for measurements in accordance with the manufacturer's instructions (Figure 2)*

**Fig. 1**

**Fig. 2**

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**EN**

- *With the nozzle of the probe approximately 2.5 cm from the surface, scan the downstream side of the HEPA filters, by passing the photometer probe in slightly overlapping strokes over the entire surface.*

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
**Fig. 3**

**EN**

- Also scan the entire periphery of the filter and the joint between filter and filter mounting frame at a scanning rate that does not exceed 5cm/s.

**Fig. 4**

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**EN CANOPY**

- In case of canopy are installed and can't remove it, scan the filter by the exhaust removing the registry cover canopy.  
**GS has not been tested for the canopy. It must be tested after installation.**

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
**EN 5. ACCEPTANCE CRITERIA**

- The photometer signal after the filter should be smaller or equal to a 0,01 % of the specific concentration before the filter.

**Figura 5.** Resultado correcto: no fuga  
*Correct result: there isn't Leak*

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**EN 1. PURPOSE**

- Check that the speed control in different points of the working area is suitable to achieve product protection. To determine the speed average as well as the uniformity range in the laminar flow zone.

**2. SCOPE**

- Bio II Advance

**3. TOOLS**

- The tool to use is an anemometer.

**4. METHODOLOGY**

- Inside the considered working area according to each cabinet model divide the working area in squares (see attached table). Then switch on the cabinet and wait until the velocity and the temperature
- Place the anemometer at 270mm high of the shelves, over each square (take all readings with the anemometer in the same direction and without covering any holes, whenever possible).
- Take note of the readings and the position of the anemometer in each point.

For this calculation, we proceed to the following formula:

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**EN 5. ACCEPTANCE CRITERIA**

- *The average of the airflow velocity at the work surface should be 0,35 m/s  $\pm$  20%.*

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- V. del anemómetro:
- V. del display:

**EN 1. PURPOSE**

- Check correct work of the alarms whwn flow speed is minor or major than the admissible.

**2. RANGE**

- Bio II Advance*

**3. TOOLS**

- Omnidirectional anemometer

**4. METHODOLOGY**

- Accede to Technical Service menu through *Password Petition* screen (User's Manual chapter. 4.5 Point 1), using for it random security password that you ask for.

- Move to *Set Points Impulsion* screen and wait some seconds until the stabilization of the cabinet.

- Choose work area's lower lecture obtained in *Downflow velocity Test* (Annex A9.3) (avoid near zones to the cabinet walls). Place the anemometer at that point and decrease slowly downflow velocity through function

button until "Insufficien Laminar Flow" alarm activated. Take note:

- Anemometer V. :
- Display V.:

- Repeat the process with higher lecture and increasing slowly downflow velocity through

the function button until "Excessive Laminar Flow" activated. Take note:

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**EN****5. ACCEPTATION CRITERIA**

- Anemometer reading at point of lowest speed  $\geq 0.28$  m/s.
- Anemometer reading at point of highest speed  $\leq 0.42$  m/s.

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**EN 1. PURPOSE**

- The purpose of this test is to calculate the inflow velocity through the front aperture of the equipment to obtain the Class II to the operator protection.

**2. SCOPE**

- Bio II Advance

**3. TOOLS**

- The tool to use is an anemometer.

**4. METHODOLOGY**

The following process shall be performed while the equipment is working:

- Take 5 velocity readings at 2,5 cm of the protective grills over the exhaust filter.
- Calculate the average velocity:  $V_{average} (m/s) = 1/5 \times (V_1 + V_2 + V_3 + V_4 + V_5)$
- Calculate the flow exhaust:  $Q_{exhaust} (m^3/h) = (V_{average} / Exhaust surface) \times 3600$

NOTE: Exhaust surface: 0,078 m<sup>2</sup>.

- Finally, calculate the inflow velocity by the following formula:

$$V_{INFLOW} (m/s) = Q_{exhaust} / (Front Aperture surface \times 3600)$$

NOTA: Inflow surfaces:

Bio II Advance 3: 0,149 m<sup>2</sup>  
 Bio II Advance 4: 0,203 m<sup>2</sup>  
 Bio II Advance 6: 0,313 m<sup>2</sup>

**EN 5. ACCEPTANCE CRITERIA**

- Test results are considered to be satisfactory when the inflow velocity is suitable to achieve operator protection ( $\geq 0,55$  m/s).

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**EN 1. PURPOSE**

- Check the correct operation of inflow air velocity alarm when flow is minor than the admisible.

**2. RANGE**

- *Bio II Advance*

**3. METHODOLOGY**

- Accede to Technical Service menu through *Password Petition* screen (User's Manual Chapter 4.5 Point 1), using for it random security password that you ask for.
- Move to *Set Points Impulsion* screen and wait some seconds until the stabilization of the cabinet
- Decrease slowly entry velocity through the

function button until "Insufficient Extraction flow" alarm is activated. Take note:

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**EN 1. OBJECT**

- *Verify the cabinet for personnel and product protection.*

**2. SCOPE**


- *Bio II Advance*

**3. TOOLS**


- *Smoke generator*

**4. METHODOLOGY**

- *Product protection: smoke is produced 5-10 cms from the protective glass (Fig. 1 and 2). The behaviour of the smoke is checked visually.*
- *Smoke is produced inside the cabinet, as close as possible to the lower edge of the glass (Fig. 3). The behaviour of the smoke is checked visually.*
- *User protection: smoke is produced inside the cabinet, as close as possible to the lower edge of the glass (Fig. 4). The behaviour of the smoke is checked visually.*

**Fig. 1**

**Fig. 2**

**Fig. 3**

**Fig. 4**

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## **EN 5. ACCEPTATION CRITERIA**

*Product protection:*

- *The generated smoke in the outside of the front aperture will be dragged to the inflow area. No fume sign have to reach the working surface.*
- *Turbulences in the working area shall not be created. The smoke generated should be dragged by the laminar flow to the surface of the working area.*

*User protection:*

- *No smoke shall be refluxed out of the cabinet once drawn in, guaranteeing the protection of the operator and environment.*

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**EN 1. OBJECT**

*Check that the level and uniformity of the lighting system is the correct one to guarantee a safety working area inside the cabinet.*

**2. SCOPE**

*Bio II Advance*

**3. TOOLS**

*Light meter.*

**4. METHODOLOGY**

- *An imaginary division of the working surface will be done in equal parts (not more than 8, neither less than 6)*
- *The illumination level will be measured in each of these parts, placing the light meter over the working area*
- *The average of the measures carried out will be calculated as well as the safety cabinet illumination level*

**Fig. 1**

Punto/ Point/	1	2	3	4	5	6	7	8	Promedio/ Average/
Valor/ Value/									

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**EN 1. ACCEPTATION CRITERIA**

*The test will be correct if:*

- *The illumination level average is **not** less than 1000 lux.*
- *Any individual reading can be less than the 70% of the average illumination level.*

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## **GUARANTEE CERTIFICATE**

*TELSTAR TECHNOLOGIES, S.L. guarantees the substitution of all the defect pieces in the period of TWO YEARS, starting in the date of the expedition of the equipment.*

*This equipment has undergone satisfactorily all the tests and revisions. Therefore, any damage caused by transport, wrong installation or inappropriate use is excluded of this guarantee.*

*This guarantee does not include any of the fungible materials comprised in the equipment, such as filters, fluids, filaments, etc.*

*The inappropriate use of the equipment excludes Telstar of any responsibility for eventual direct or indirect damages of any nature.*

*Likewise, this guarantee will be void if pieces other than the ones supplied by TELSTAR TECHNOLOGIES, S.L. are used for the repairing of the equipment.*

**NOTE. This guarantee is valid for all the equipment unless otherwise stated. In order to obtain any information about this guarantee, please contact with our Sales Department.**

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