


SensoLyt[®] 900-P



MultiLine[®] INTELLIGENT DIGITAL SENSORS

IDS pH electrode with polymer electrolyte

Operating manual

Copyright

© Weilheim 2012, WTW GmbH
Reprinting - even as excerpts - is only allowed with the explicit written
authorization of WTW GmbH, Weilheim.
Printed in Germany.

8 General information

Automatic sensor recognition

The sensor electronics with the stored sensor data are in the connecting head of the electrode. The data include, among other things, the sensor type and series number. In addition, the calibration data are stored in the sensor with each calibration and the calibration history is recorded (the last 10 calibrations). The data is recalled by the meter when the sensor is connected and is used for measurement and for measured value documentation.

Storing the calibration data in the sensor ensures that the correct slope and asymmetry are automatically used if the sensor is operated with different meters. On the other hand, different calibrated sensors can be used with one meter without the need to recalibrate.


The digital transmission technique guarantees the failure-free communication with the meter even with long connection cables. If the sensor firmware is enhanced by WTW, it can be updated via the meter.

Armoring A 925-P/K and A 925-P/S (accessories)

To protect the electrode against mechanical impacts, the A 925-P/K armoring (with plastic protective hood) or A 925-P/S (with steel protective hood) can be installed and uninstalled as necessary. The steel model is simultaneously used as a sinker for depth measurements.

9 Commissioning, measuring, calibration

9.1 Opening and closing the IDS plug connection



Opening the plug connection

- If necessary, clean the plug connection
- Open the locking device (step 1)
- Use your thumb and index finger to press the clips of the connector together, and pull the connector out of the plug (step 2).

Closing the plug connection

- Make sure that the plug connection is completely dry and clean.
- Align the guidance of the connector with the groove in the plug and insert the connector in the unlocked plug until it catches (step 3).
- Close the locking device (step 4).

9.2 Commissioning

Prepare the electrode for measuring as follows:

- Remove the watering cap from the electrode tip. Possible salt deposits in the area of the watering cap do not affect the measuring characteristics and can easily be removed with deionized water.



Please keep the watering cap. It is required for the electrode to be stored. Always keep the watering cap clean.

- Connect the electrode to a free IDS sensor plug-in position of the multi parameter probe or to an IDS connection of the meter.
To open and close the IDS plug-in position please note the section 9.1 OPENING AND CLOSING THE IDS PLUG CONNECTION.




Connection cables in different lengths to connect the electrode to the meter are listed in section 13 WEAR PARTS AND ACCESSORIES.

- Calibrate the electrode according to the operating manual of the meter and observe the rules in the following section while doing so:

9.3 Calibration and measurement: General rules

- Avoid the carryover of any solution (sample or buffer solution) from one measurement to the next by taking the following measures:
 - Shortly rinse the calibration and sample beakers with the solution the beakers are to be filled with next.
 - Between measurements, rinse the electrode with the solution that follows. Alternatively, you can also rinse the electrode with deionized water and then carefully dab it dry.
- To measure in aqueous solutions, it is recommended to immerse the electrode in a vertical or slightly tilted position.
- Observe the correct depth of immersion and make sure the contact between the junction and test sample is thorough. The junction (hole junction) is in the area of the bottom end of the shaft.



- For measurements in aqueous solutions, provide approximately the same stirring conditions for measuring as for calibrating.

NOTE

Prevent contact of the pH membrane to the beaker bottom to avoid scratches on the pH membrane.

Subsequent calibrations

The frequency of subsequent calibrations depends on the application. The meters provide an option where you can enter a calibration interval. After the calibration interval has expired, the meter will automatically remind you of the due calibration.

10 Storage

During short measuring breaks

Immerse the electrode in reference electrolyte (KCl 3 mol/L, Ag⁺ free). Prior to the next measurement, shortly rinse the electrode with the test sample or deionized water.

Overnight or longer

Put the clean electrode in the watering cap that is filled with reference electrolyte (KCl 3 mol/L, Ag⁺ free).

NOTE

Do not store the electrode dry or in deionized water.

The electrode could be permanently damaged by this.

If the liquid in the watering cap has dried up, condition the electrode in reference electrolyte (KCl 3 mol/L, Ag⁺ free) for at least 24 hours.



During longer storing periods, salt sediments may develop on the watering cap. They do not affect the measuring characteristics and can easily be removed with deionized water when the electrode is put into operation again.

11 Aging

Every pH electrode undergoes a natural aging process. With aging, the responding behavior becomes slower and the electrode slope and asymmetry change. Moreover, extreme operating conditions can considerably shorten the lifetime of the electrode. These are:

- Strong acids or lyes, hydrofluoric acid, organic solvents, oils, fats, bromides, sulfides, iodides, proteins
- High temperatures
- High changes in pH and temperature.

The warranty does not cover failure caused by measuring conditions and mechanical damage.

12 Maintenance and cleaning

Cleaning

Remove water-soluble contamination by rinsing with deionized water. Other types of contamination have to be removed as follows while the contact time with the detergents should be kept as short as possible:

Contamination	Cleaning procedure
Fat and oil	Rinse with water containing household washing-up liquid
Lime and hydroxide deposits	Rinse with citric acid (10 % by weight)

NOTE

Hydrofluoric acid, hot phosphoric acid and strong alkaline solutions destroy the glass membrane.

After cleaning

Rinse the electrode with deionized water and condition it in reference electrolyte solution for at least 1 hour. Then recalibrate the electrode.

13 Wear parts and accessories

Maintenance equipment

Description	Model	Order no.
Reference electrolyte solution 250 ml to fill the watering cap (KCl 3 mol/l, Ag ⁺ -free)	KCl-250	109 705

Connection cable SensoLyt® 900-P - meter

Description	Model	Order no.
IDS connection cable, 1.5 m	AS/IDS-1.5	903 850
IDS connection cable, 3 m	AS/IDS-3	903 851
IDS connection cable, 6 m	AS/IDS-6	903 852
IDS connection cable, 10 m	AS/IDS-10	903 853
IDS connection cable, 15 m	AS/IDS-15	903 854
IDS connection cable, 20 m	AS/IDS-20	903 855
IDS connection cable, 25 m	AS/IDS-25	903 856
IDS connection cable, 40 m	AS/IDS-40	903 857
IDS connection cable, 60 m	AS/IDS-60	903 858
IDS connection cable, 100 m	AS/IDS-100	903 859

General accessories

Description	Model	Order no.
Blind plug for IDS plug	BPO/IDS 900	908 371
Armoring without protective hood	A 925-P	903 838
Armoring with plastic protective hood	A 925-P/K	903 839
Armoring with steel protective hood	A 925-P/S	903 840



Detailed information on our wide range of buffer solutions and more accessories is given in the price list of the WTW catalog "Laboratory and field instrumentation".

14 Technical data

General features	Reference electrolyte	Polymer electrolyte, AgCl free
	Junction	1-hole junction
	Shunt conduction element	Ag/AgCl
	Temperature sensor	Integrated NTC 30 (30 kΩ at 25 °C / 77 °F)
Measurement and application characteristics	pH measuring range	0.000 ... 12.000
	Allowed temperature range	0 ... 60 °C
	Membrane resistance at 25 °C	< 500 MOhm
	Typical application	Field
Accuracy of the IDS measuring technique	Measured parameter	Accuracy (± 1 digit)
	pH	± 0.004
	U [mV]	± 0.2
	T [°C]	± 0.1
Pressure range at temperature	<u>Temperature</u>	<u>allowed overpressure</u>
	0 °C (32 °F)	1000 kPa (10 bar)
	20 °C (68 °F)	1000 kPa (10 bar)
	30 °C (86 °F)	500 kPa (5 bar)
	40 °C (104 °F)	300 kPa (3 bar)
	60 °C (140 °F)	100 kPa (1 bar)
	The electrodes meet the requirements according to article 3(3) of the directive, 97/23/EC ("pressure equipment directive").	
IDS plug	Connection type	4-pole, watertight plug connection with lock, reverse polarity protected

**Shaft dimensions,
material**

Shaft length	120 mm
Shaft diameter	12 mm
Shaft material	Glass
IDS plug	<ul style="list-style-type: none">● Synthetic materials: Glass fiber reinforced Noryl, TPU, TPC-ET, POM, PVC, PEEK, PBT● O-ring: FPM● Contacts gold-plated

Storage With watering cap; filled with KCl 3 mol/L, Ag⁺ free

Disposal Residual waste



Wissenschaftlich-Technische Werkstätten GmbH

Dr.-Karl-Slevogt-Straße 1
D-82362 Weilheim

Germany

Tel: +49 (0) 881 183-0
+49 (0) 881 183-100
Fax: +49 (0) 881 183-420
E-Mail: Info@WTW.com
Internet: <http://www.WTW.com>