**Instruction Manual**  
**Manuale di istruzioni**  
**Manuel d'instructions**  
**Manual de instrucciones**  
**Bedienungsanleitung**  
**指导手册**

## RESPIROMETRIC Sensor

F102B0153

### General Information / Informazioni Generali / Informations Générales / Información General / Allgemeine Hinweise / 一般信息

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#### **This unit has been designed and manufactured in compliance with the following standards:**

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Safety requirements for electrical equipment for measurement, control and for laboratory use

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We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following standards:

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Machinery directive 2006/42/EC / Macchine 2006/42/CE / Machines 2006/42/CE / Máquinas 2006/42/CE /

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## Declaration of conformity

We, the manufacturer VELP Scientifica, under our responsibility declare that the product is manufactured in conformity with the following regulations:

S.I. 2016/1101 Electrical Equipment (Safety) Regulations 2016

S.I. 2016/1091 Electromagnetic Compatibility Regulations 2016

according to the relevant designated standards:

EN 61010-1 Safety requirements for electrical equipment for measurement, control, and laboratory use  
Part 1: General requirements

EN 61326-1 Electrical equipment for measurement, control and laboratory use - EMC requirements -  
Part 1: General requirements

and satisfies the essential requirements of regulations:

S.I. 2008/1597 Supply of Machinery (Safety) Regulations 2008

S.I. 2012/3032 Restriction of the Use of Certain Hazardous Substances in Electrical and Electronic Equipment  
Regulations 2012

S.I. 2013/3113 Waste Electrical and Electronic Equipment Regulations 2013

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## 1. Introduction

BOD measurement with RESPIROMETRIC Sensor is based on the RESPIROMETRIC method that provides a direct measurement of the oxygen consumed by microorganisms from an air or oxygen-enriched environment in a closed atmosphere, at a constant temperature ( $20\text{ °C} \pm 0.5$ ) and with a soft sample stirring.

The measurement of BMP (Biochemical Methane Potential or Biochemical Methanization Potential) represents the anaerobic correspondent of the BOD tests in fact an adequate anaerobic microbial inoculation is used to degrade the test substrate anaerobically. During the test the production of biogas and/or methane is proportional to the quantity of organic substance gradually degraded.

## 2. Assembly and installation

Upon receipt and after having removed the packaging, please check the integrity of the instrument. The box includes:

- RESPIROMETRIC Sensor
- Alkali holder
- Instruction manual

The RESPIROMETRIC Sensor is delivered with 2 lithium batteries installed: remove the insulating strip from the rear without opening it.

## 3. Operating controls

The display of the electronic measuring unit is normally turned off reducing the consumption of internal batteries. The touch of a key activates the display which automatically is turned off after three seconds. The key activates different functions if the unit is performing a measurement or if it is in memory code and if the unit is in wireless mode. When the display is turned off and no analysis is performed, RESPIROMETRIC Sensor shows the symbol  $\square\square$  by pressing any keys, in wireless mode.

### 3.1 Wireless or non-wireless mode

RESPIROMETRIC Sensor is already set on the wireless mode. To switch to a non-wireless mode, keep pressed SET when the display is off until "S.E.T." appears. When the three dots disappear:

- Key **SET** Pass to wireless mode ON
- Key **START** Pass to wireless mode OFF

### 3.2 Description of keys (wireless mode)

Out of a measuring cycle (before starting)

- Key **SET** with display off Wireless symbol  $\square\square$  is shown
- Key **SET** with display on Wireless symbol  $\square\square$  is shown
- Key **START** with display off Wireless symbol  $\square\square$  is shown
- Key **START** with display on The cycle of measurement is started

During a measuring cycle (after starting)

- Key **SET** The scale in use is shown
- Key **START** The current measure is shown (also after the 5-day period if Reset is not performed)
- Key **START** kept pressed The memory mode is entered (also after the 5-day period if Reset is not performed)

Memory mode (it is active during a measuring cycle and also after if Reset is not performed)

### 3.3 Description of keys (non-wireless mode)

Out of a measuring cycle (before starting)

- Key **SET** with display off The scale is shown
- Key **SET** with display on The scale is changed
- Key **START** with display off The scale is shown
- Key **START** with display on The cycle of measurement is started

During a measuring cycle (after starting)

- Key **SET** The scale is shown
- Key **START** The current measure is shown (also after the 5-day period if Reset is not performed)
- Key **START** kept pressed The memory mode is entered (also after the 5-day period if Reset is not performed)

Memory mode (it is active during a measuring cycle and also after if Reset is not performed)

- Key **SET** The stored value is shown
- Key **START** The order number of stored value is increased

If no values are stored, the display shows three points after a prolonged pressing of key START.

Reset

Keys **SET+START** kept pressed Measurement cycle stopped, stored values deleted and the scale is shown. Before starting a new measuring cycle, it is requested to reset the stored values.

**NOTE:** before beginning an additional RESPIROMETRIC Sensor, it's necessary to wait that the display shows "An".

### 3.4 Last data recovery

If during an analysis the Wireless DataBox™ loses some data, it's possible to retrieve the data for the last seven samplings by pressing SET twice: in this way data will be automatically sent to the Wireless DataBox™.

### 3.5 Scale selection

Scale	Volume of sample
0 ÷ 4000 mg O <sub>2</sub> /l	25ml
0 ÷ 999 mg O <sub>2</sub> /l	100 ml
0 ÷ 600 mg O <sub>2</sub> /l	150 ml
0 ÷ 250 mg O <sub>2</sub> /l	250 ml
0 ÷ 90 mg O <sub>2</sub> /l	400 ml
PRE	Refer to the method
BMP	Refer to the method

The scale is related to the expected value. The possibility of measuring diluted samples extends the range. Results are shown as mg O<sub>2</sub>/l (ppm) for BOD test and values in mbar for PRE scale.

#### Wireless mode

Before the analysis set the scale on the RESPIROSoft™

#### Non-wireless mode (after the reset)

Key **SET** (display on)                      The scale is changed  
Key **SET** or **START**                         The scale is shown

### 3.6 Analysis duration and sampling time

#### Wireless mode

With RESPIROMETRIC Sensor an analysis can last from 5 to 120 days. The sampling time can be selected depending on the desired duration of the analysis, from 1 h to 2 days

#### Non-wireless mode

Respirometric Sensor will perform only 5-day test with 24-hour sampling time for BOD analysis or 30 days using PRE scale.

#### Wireless mode

Before the analysis set the scale on the RESPIROSoft™

### 3.7 Alarms

- AL1** Calibration error. Restart the analysis. If the problem persists contact the service.
- AL2** Analysis start denied. Check on RESPIROSoft™.
- AL3** Max. sensors q.ty reached. 50 analyses are already in progress. Wait the end.
- AL4** One parameter or more are wrong
- AL5** No PC response. Check the distance between the RESPIROMETRIC Sensor and the Wireless DataBox™.
- AL6** The analysis is stopped from the PC
- AL7** Unknown command result
- AL8** Data sending failed. Restart the analysis. If the problem persists contact the service.
- AL9** No PC response. Check if the RESPIROSoft™ is active.
- AL10** No automatic measure is done
- AL11** Wireless DataBox™ not found. Check if the Wireless DataBox™ is connected.
- AL12** Wireless DataBox™ connection isn't strong enough
- AL13** Analysis start is denied due to low battery level
- bAt** Low battery (wireless mode OFF). Remove the batteries with the fingers and replace them with new ones paying attention to the polarity indicated inside the instrument.
- Err** Calibration error (wireless mode OFF).
- The measured value overcomes the maximum value of the chosen scale.

## 4. Maintenance

RESPIROMETRIC Sensor warns about the required change of batteries when the unit is resetted and the charge is lower than a minimum threshold. Use a coin lever in the rear part of RESPIROMETRIC Sensor and open it. Batteries are removed using the fingers. Introduce two fresh batteries at the same time paying attention to the polarity. Battery replacement deletes any stored value.

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A00000136

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## Thank you for having chosen VELP!

Established in 1983, VELP is today one of the world's leading manufacturers of analytical instruments and laboratory equipment that has made an impact on the world-wide market with Italian products renowned for innovation, design and premium connectivity. VELP works according to **ISO 9001**, **ISO14001** and **ISO 45001** Quality System Certification.

Our instruments are manufactured in Italy according to the IEC 1010-1 and CE regulation.

Our product lines:

### Analytical instruments

Elemental Analyzers  
Digestion Units  
Distillation Units  
Solvent Extractors  
Fiber Analyzers  
Dietary Fiber Analyzers  
Oxidation Stability Reactor  
Consumables

### Laboratory Equipment

Magnetic Stirrers  
Heating Magnetic Stirrers  
Heating Plates  
Overhead stirrers  
Vortex mixers  
Homogenizers  
COD Thermoreactors  
BOD and Respirometers  
Cooled Incubators  
Flocculators  
Overhead Shakers  
Turbidimeter  
Radiation Detector  
Open Circulating Baths  
Pumps

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