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Tootja





REF Referents = toote number

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Operation Instructions: 'haematocrits'

9µL, 18µL, 60µL, 70µL, 75µL, 100µL

(Sodium heparin, ammonium heparin, without coating)

Instructions for use

The capillaries are disposable glass capillary pipettes intended for one-time use with subsequent disposal.

Suitable for medical professionals only.

Application temperature: +2 – +40°C

National guidelines for occupational safety and quality assurance must be observed.

Purpose

Intended purpose: Haematocrits for transferring defined heparin of capillary whole blood for quantitative analysis. The heparin contained in the glass capillary prevents the coagulation of blood in order to avoid a negative impact on the downstream analysis. The blood is entered into the capillary directly at the analysis device.

Prohibitions of use: Haematocrits are NOT intended for the following applications: Sampling of other liquids than capillary whole blood, storage of samples, measurements within the capillary itself.

Hazards and safety

The precautionary measures required for the use of laboratory reagents and body fluids as well as microbiological samples must be observed. Handling may only be carried out by qualified personnel. National and internal laboratory guidelines for occupational safety and infection control must be observed. Suitable protective clothing, protective goggles and disposable gloves must be worn during work.

Attention must be paid to effective protection against infection in accordance with the laboratory guidelines.

The user is committed to report any serious incidents in connection with the product to the competent authority of the Member State in which the user and/or the patient is registered and to Hirschmann Laborgeräte GmbH & Co. KG.

The capillaries comply with the requirements of Regulation (EU) 2017/746 regarding in vitro diagnostic medical devices.

Disposal

The national statutory regulations must always be observed when disposing of the capillaries. Used capillaries contaminated with potentially infectious substances of human origin may pose a risk of injuries, infections or microbiological hazards. Potentially contaminated capillaries must be disposed of in such a way that injuries, infections or microbiological hazards to third parties are excluded. Unused capillaries can be disposed of with the household waste. Unused capillaries must be disposed of in such a way that injuries to third parties are excluded.

No.	Description
9100170	HAEMATOCRITS OD. 1.5-1.6 MM 70 µL
910017048	HAEMATOCRITS OD. 1.5-1.6 MM 70 µL
9100160	HAEMATOCRIT CAPILLARIES 60µl OD. 1.3-1.4 MM
9100175	HAEMATOCRIT CAPILLARIES 75µl OD. 1.5-1.6 MM
910017547	HAEMATOCRITS OD. 1.5-1.6 MM 75 µL
910017548	HAEMATOCRITS OD. 1.5-1.6 MM 75 µL
910021801	HAEMATOCRIT; NA-HEP. 18 µL OD1.35 +/- 0,05; ID0,55 +/-0,05
9100260	HAEMATOCRIT CAPILLARIES 60 µl OD. 1.3-1.4 MM; NA-HEP.
9100270	HAEMATOCRITS OD. 1.5-1.6 MM 70 µL; SODIUM-HEP.
910027048	HAEMATOCRITS OD. 1.5-1.6 MM 70 µL; SODIUM HEP.
9100275	HAEMATOCRIT CAPILLARIES 75µl OD. 1v5-1.6 MM; NA-HEP.
910027547	HAEMATOCRITS OD. 1.5-1.6 MM 75 µL; SODIUM-HEP.
910027548	HAEMATOCRITS OD. 1.5-1.6 MM 75 µL; SODIUM HEP.
9100299	HAEMATOCRIT CAPILLARY 100µl OD. 1.5-1.6 MM; NA-HEP
9100418	HAEMATOCRIT CAPILLARY 18 µl OD. 1.3-1.4 MM; AMMONIUM HEP.
9100475	HAEMATOCRITS OD. 1.5-1.6 MM 75 µL; AMMONIUM HEP.
9030208	CAPILLARIES NA-HEP; 9 µl IN VIAL; LENGTH 32 MM

IVD In vitro diagnostic medical devices

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