Glucose OF Medium ISO 21528

Semisolid medium for determining oxidative or fermentative metabolism of bacteria.

DESCRIPTION

Glucose OF medium is a semisolid medium used for the identification of gram-negative bacteria on the basis of their ability to oxidize or ferment glucose.

This medium complies with ISO 21528:2017 for confirmation of Enterobacteriaceae from food, animal feeding stuffs and environmental samples.

TYPICAL FORMULA

	(g/l)
Enzymatic Digest of Casein	2.0
Glucose	10.0
Sodium Chloride	5.0
Dipotassium Hydrogen Phosphate	0.3
Bromthymol Blue	0.08
Agar	3.0

Final pH 6.8 ± 0.2 at 25°C

METHOD PRINCIPLE

Enzymatic digest of casein provides amino acids, nitrogen, carbon, vitamins and minerals for organisms growth. Glucose is the fermentable carbohydrate. Sodium chloride maintains the osmotic balance of the medium. Dipotassium phosphate acts as buffer agent. Bromthymol blue is the pH indicator. Agar is the solidifying agent.

PREPARATION

Dehydrated medium Suspend 20.38 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. Dispense into tubes. Sterilize in autoclave at 121°C for 15 minutes.

TEST PROCEDURE

Just before use, heat the medium in boiling water or flowing steam for 15 min to remove oxygen, then cool rapidly to the incubation temperature.

Using a wire, inoculate tubes by stabbing colonies into the medium. Overlay with minimal 1 cm of sterile mineral oil (Vaseline Oil Droppers, ref. 87006). Incubate at 37°C for 24 hours.

INTERPRETING RESULTS

If a yellow color develops throughout the medium the fermentation test is positive.

According to ISO 21528, colonies that are oxidase-negative and glucose-positive are confirmed as Enterobacteriaceae.

APPEARANCE

Dehydrated medium: free-flowing, homogeneous, beige to greenish.

Prepared medium: slightly opalescent, green.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store tubes at 10-25°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 4 years.

Medium in tubes: 2 years.

QUALITY CONTROL

The medium is inoculated with the microbial strains indicated in the QC table.

Incubation conditions: 37 ± 1°C for 24 ± 2 hours.

QC Table.

Microorganism		Growth	Fermentation test
<i>Escherichia coli</i>	WDCM 00012	Good	Positive (+), yellow color
<i>Pseudomonas aeruginosa</i>	WDCM 00025	Good	Negative (-), yellow color only at the top of the tube

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for professional use only and must be used by properly trained operators.

DISPOSAL OF WASTE

Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. ISO 21528-1:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of Enterobacteriaceae – Part 1: Detection of Enterobacteriaceae.
2. ISO 21528-2:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of Enterobacteriaceae – Part 2: Colony-count technique.
3. EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
4. Welch, D.F., et al (1987) Selective and differential medium for recovery of *Pseudomonas cepacia* from the respiratory tracts of patients with cystic fibrosis. J. Clin. Microbiol., 25; 1730-1734.
5. Costin, I.D. (1967) An outline for the biochemical identification of aerobic and facultatively anaerobic gram-negative rods of medical interest. Intern. Congr. f. Chemotherapie Wien, B2/1; 73-76.
6. Mossel, D.A.A. Martin, G. (1961) Milieu simplifié permettant l'étude des divers modes d'action des bactéries sur les hydrates des carbone. - Ann. Inst. Pasteur de Lille, 12; 225-226.
7. Hugh, R., A. Leifson, E. (1953) The taxonomic significance of fermentative versus oxidative metabolism of carbohydrates by various gram-negative bacteria. J. Bact., 66; 24-26.

PRESENTATION

		Contents	Ref.
Glucose OF Medium ISO 21528	Tubes	20 x 10 ml tubes	24482
Glucose OF Medium ISO 21528	Dehydrated medium	500 g of powder	610388
Glucose OF Medium ISO 21528	Dehydrated medium	100 g of powder	620388
Glucose OF Medium ISO 21528	Dehydrated medium	5 kg of powder	6103885

TABLE OF SYMBOLS

LOT Batch code	 Do not reuse	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	

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BIBLIOGRAFIA

1. ISO 21528-1:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of Enterobacteriaceae – Part 1: Detection of Enterobacteriaceae.
2. ISO 21528-2:2017. Microbiology of the food chain – Horizontal method for the detection and enumeration of Enterobacteriaceae – Part 2: Colony-count technique.
3. EN ISO 11133:2014. Microbiology of food, animal feed and water – Preparation, production, storage and performance testing of culture media.
4. Welch, D.F., et al (1987) Selective and differential medium for recovery of *Pseudomonas cepacia* from the respiratory tracts of patients with cystic fibrosis. J. Clin. Microbiol., 25; 1730-1734.
5. Costin, I.D. (1967) An outline for the biochemical identification of aerobic and facultatively anaerobic gram-negative rods of medical interest. Intern. Congr. f. Chemotherapie Wien, B2/1; 73-76.
6. Mossel, D.A.A. Martin, G. (1961) Milieu simplifié permettant l'étude des divers modes d'action des bactéries sur les hydrates de carbone. - Ann. Inst. Pasteur de Lille, 12; 225-226.
7. Hugh, R., A. Leifson, E. (1953) The taxonomic significance of fermentative versus oxidative metabolism of carbohydrates by various gram-negative bacteria. J. Bact., 66; 24-26.

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		Contenuto	Ref.
Glucose OF Medium ISO 21528	Provette	Provette 20 x 10 ml	24482
Glucose OF Medium ISO 21528	Terreno disidratato	500 g di polvere	610388
Glucose OF Medium ISO 21528	Terreno disidratato	100 g di polvere	620388
Glucose OF Medium ISO 21528	Terreno disidratato	5 k g di polvere	6103885

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