DNASE Agar (Dehydrated)

REF CM0321B

EN

Intended Use

DNASE Agar (Dehydrated) (CM0321B) is a medium used for the differentiation of microorganisms from clinical samples based on deoxyribonuclease activity, particularly *Staphylococcus* species.

The device is for professional use only, is not automated, nor is it a companion diagnostic.

Summary and Explanation

Staphylococci are Gram-positive cocci, widespread in nature and generally non-pathogenic^{1,2}. However, some staphylococci can cause disease in humans^{1,2}, most frequently *S. aureus*², an opportunistic pathogen commonly found as a commensal of the human upper respiratory tract and skin microbiomes¹. *S. aureus* enters the body through cuts or abrasions where it can cause a spectrum of disease, from minor infections such as cellulitis and abscesses, to life-threatening illnesses such as pneumonia, toxic shock syndrome and sepsis¹.

S. aureus is one of the leading causes of hospital acquired infections, further complicated by the rapid emergence of antibiotic resistance in this pathogen. Strains of *S. aureus* that have developed multiple drug-resistance to beta-lactam antibiotics, termed methicillin-resistant *S. aureus* (MRSA), are a significant problem in hospitals as well as the community³.

There are two *Serratia* species which are commonly isolated from clinical samples, including *Serratia marcescens*. Members of this genus characteristically produce three enzymes; lipase, DNase and gelatinase. *Serratia* species are found in faeces, wound exudates, respiratory specimens, blood, eye cultures, and urine⁴.

Principle of Method

Tryptose supplies nitrogen and other growth factors to support the growth of bacteria. Deoxyribonucleic acid is the substrate for DNase, sodium chloride is added to maintain the osmotic equilibrium and agar is the solidifying agent.

Typical Formula

	<u>grams per litre</u>
Tryptone	20.0
Deoxyribonucleic acid	2.0
Sodium chloride	5.0
Agar	12.0

Materials Provided

500g of dehydrated DNase Agar that yields approximately 12.8L after reconstitution.

Materials Required but Not Supplied

- Inoculating loops, Incubators
- Quality control organisms
- Sterile containers
- 1N hydrochloric acid

Storage

- Store product in its original packaging between 10°C and 30°C.
- Keep container tightly closed.
- The product may be used until the expiry date stated on the label.
- Protect from moisture.
- Store away from light.
- Allow reconstituted product to equilibrate to room temperature before use.

Once reconstituted, store media between 2°C and 10°C.

Warnings and Precautions

- For in vitro diagnostic use only.
- For professional use only.
- Inspect the product packaging before first use.
- Do not use the product if there is any visible damage to the packaging (pot or cap).
- Do not use the product beyond the stated expiry date.
- Do not use the device if signs of contamination are present.
- It is the responsibility of each laboratory to manage waste produced according to their nature and degree of hazard and to have them treated or disposed of in accordance with any federal, state and local applicable regulations. Directions should be read and followed carefully. This includes the disposal of used or unused reagents as well as any other contaminated disposable material following procedures for infectious or potentially infectious products.

- Ensure the lid of the container is kept tightly closed after first opening and between use to minimise moisture ingress, which may result in incorrect product performance.

Refer to the Safety Data Sheet (SDS) for safe handling and disposal of the product (www.thermofisher.com).

Serious Incidents

Any serious incident that has occurred in relation to the device shall be reported to the manufacturer and the relevant regulatory authority in which the user and/or the patient is established.

Specimen Collection, Handling and Storage

Specimens should be collected and handled following local recommended guidelines, such as the UK Standards for Microbiology Investigations (UK SMI) ID 7, ID 16, TP 12, B 29 and Q 5.

Procedure

Suspend 39g in 1 litre of distilled water. Bring to the boil to dissolve completely. Sterilise by autoclaving at 121°C for 15 minutes. Cool to 50°C. Mix well and pour into sterile containers.

Interpretation

Once the medium is flooded with acid, *Staphylococcus aureus* and *Serratia marcescens* will appear as growth with clear zones. *Staphylococcus epidermidis* and *Klebsiella pneumoniae* will appear as growth with no clearing.

Quality Control

It is the responsibility of the user to perform quality control testing taking into account the intended use of the medium, and in accordance with any local applicable regulations (frequency, number of strains, incubation temperature etc.).

The performance of this medium can be verified by testing the following reference strains.

Incubation Conditions: 18h at 37°C aerobic

Positive Controls	
Inoculum level: 10 ⁴ – 10 ⁶ cfu	
<i>Staphylococcus aureus</i> ATCC® 25923™	Growth with clear zones
<i>Staphylococcus aureus</i> ATCC® 9144™	Growth with clear zones
<i>Serratia marcescens</i> ATCC® 8100™	Growth with clear zones
<i>Staphylococcus epidermidis</i> ATCC® 12228™	Growth, no clearing
<i>Staphylococcus epidermidis</i> ATCC® 14990™	Growth, no clearing
<i>Klebsiella pneumoniae</i> ATCC® 13883™	Growth, no clearing

Limitations

The DNase reaction for staphylococci is an indication of pathogenicity, it cannot be used as the sole criterion for identification; additional tests must be used. Both *S. aureus* and *S. epidermidis* produce DNase, however the former produces it in larger quantities. *Staphylococcus intermedius* and *schleiferi* may also produce DNase. Small zones of clearing may be caused by other enzymes or organic acid production. Organisms other than staphylococci, *serratia* and *aeromonads* can produce DNases. Once the hydrochloric acid has been applied to the medium the plate must be read within a few minutes and further testing cannot be carried out by reincubation. Normal hydrochloric acid is bactericidal and test organisms cannot be recovered from the surface of the agar after flooding.

Performance Characteristics

Accuracy has been demonstrated through review of the QC data. Correct detection of DNase positive strains of *Staphylococcus* spp. is confirmed by the inclusion of well-characterised isolates in the QC processes performed as part of the manufacture of each batch of the device. The precision of DNASE Agar (Dehydrated) (CM0321B) device was demonstrated by an overall pass rate of 100% obtained for the product over many years of testing (1995 and 2022; 306 batches of CM0321B). This shows that the performance is reproducible.

The DNASE Agar (Dehydrated) (CM0321B) device has been tested in-house as part of the QC process since the product was launched in 1995. For target organisms, when using an inoculum of 10⁴ – 10⁶ cfu of *Staphylococcus aureus* ATCC® 25923™, *Staphylococcus aureus* ATCC® 9144™ and *Serratia marcescens* ATCC® 8100™, incubating the device at 37 °C for 18 hours and flooding the agar surface with 1N hydrochloric acid, the user can effectively identify zones of clearing around deoxyribonuclease positive organisms as listed in this document. For deoxyribonuclease negative organisms, when using an inoculum of 10,000-100,000 cfu *Staphylococcus epidermidis* ATCC® 12228™, *Staphylococcus epidermidis* ATCC® 14990™

and *Klebsiella pneumoniae* ATCC® 13883™, incubating the device at 37 °C for 18 hours and flooding the agar surface with 1N hydrochloric acid, the user can effectively identify no zones of clearing around the growth, as listed in this document.

Bibliography

1. Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston.
<https://www.ncbi.nlm.nih.gov/books/NBK8448/>.
2. Public Health England. 2020a. 'Investigation of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
3. Public Health England, 2020b. 'Investigation of specimens for screening for MRSA'. UK Standards for Microbiology Investigations. B 29 Issue 7.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887595/B_29i7.pdf
4. Public Health England, 2015. 'Identification of Enterobacteriaceae'. UK Standards for Microbiology Investigations. ID 16 Issue 4.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf

Symbol Legend

Symbol	Definition
	Catalogue number
	In Vitro Diagnostic Medical Device
	Batch code
	Temperature limit
	Use-by date
	Keep away from sunlight
	Do not re-use
	Consult instructions for use or consult electronic instructions for use
	Contains sufficient for <n> tests
	Do not use if packaging damaged and consult instructions for use
	Manufacturer
	Authorized representative in the European Community/ European Union
	European Conformity Assessment
	UK Conformity Assessment
	Unique device identifier

	<p>Importer - To indicate the entity importing the medical device into the locale. Applicable to the European Union</p>
<p>Made in the United Kingdom</p>	<p>Made in the United Kingdom</p>

©2022 Thermo Fisher Scientific Inc. All rights reserved.
 ATCC and ATCC catalogue marks are a trademark of American Type Culture Collection.
 All other trademarks are the property of Thermo Fisher Scientific Inc. and its subsidiaries.

Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, UK

For technical assistance please contact your local distributor.

Revision information

Version	Modifications Introduced
2.0	2024-06-14

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Litteratur

1. Foster, Timothy J. 1996. "Staphylococcus." I Medical Microbiology. 4th Edition., 4th ed. Galveston.
<https://www.ncbi.nlm.nih.gov/books/NBK8448/>.
2. Public Health England. 2020a. 'Investigation of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 udgave
4. https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07i4.pdf
3. Public Health England, 2020b. 'Investigation of specimens for screening for MRSA.' UK Standards for Microbiology Investigations. B 29 Issue 7.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887595/B_29i7.pdf
4. Public Health England, 2015. 'Identification of Enterobacteriaceae.' UK Standards for Microbiology Investigations. ID 16 Issue 4.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/TP_16i4.pdf

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1. Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston. <https://www.ncbi.nlm.nih.gov/books/NBK8448/>.
2. Public Health England. 2020a. 'Investigation of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
3. Public Health England, 2020b. 'Investigation of specimens for screening for MRSA'. UK Standards for Microbiology Investigations. B 29 Issue 7. https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887595/B_29i7.pdf
4. Public Health England, 2015. 'Identification of Enterobacteriaceae'. UK Standards for Microbiology Investigations. ID 16 Issue 4. https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf

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Oxoid Limited, Wade Road, Basingstoke, RG24 8PW, H.B.

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Kirjandus

1. Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston.
<https://www.ncbi.nlm.nih.gov/books/NBK8448/>.
2. Public Health England. 2020 a. 'Investigation of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
3. Public Health England, 2020b. 'Investigation of specimens for screening for MRSA'. UK Standards for Microbiology Investigations. B 29 Issue 7.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887595/B_29i7.pdf
4. Public Health England, 2015. 'Identification of Enterobacteriaceae'. UK Standards for Microbiology Investigations. ID 16 Issue 4.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf

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Bibliografi

1. Foster, Timothy J. 1996. "Staphylococcus." I Medical Microbiology. 4. utgave, 4. utg. Galveston.
<https://www.ncbi.nlm.nih.gov/books/NBK8448/>.
2. Public Health England. 2020a. "Investigation of Staphylococcus species, Micrococcus species and Rothia species." UK Standards for Microbiology Investigations. ID 7 Utgave 4.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
3. Public Health England, 2020b. "Investigation of specimens for screening for MRSA." UK Standards for Microbiology Investigations. B 29 utgave 7.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887595/B_29i7.pdf
4. Public Health England, 2015. "Identification of Enterobacteriaceae." UK Standards for Microbiology Investigations. ID 16 utgave 4.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf

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Bibliografia

1. Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston. (4ª Edição, 4ª Ed. Galveston.) <https://www.ncbi.nlm.nih.gov/books/NBK8448/>.
2. Public Health England. 2020a. 'Investigation of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4. (Edição ID 7 4.) https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
3. Public Health England, 2020b. 'Investigation of specimens for screening for MRSA'. UK Standards for Microbiology Investigations. B 29 Issue 7. (Edição B 29 7.) https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887595/B_29i7.pdf
4. Public Health England, 2015. 'Identification of Enterobacteriaceae'. UK Standards for Microbiology Investigations. ID 16 Issue 4. (Edição ID 16 4.) https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf

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
Zdroje

1. Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston.
<https://www.ncbi.nlm.nih.gov/books/NBK8448/>
2. Public Health England. 2020a. 'Investigation of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
3. Public Health England, 2020b. 'Investigation of specimens for screening for MRSA'. UK Standards for Microbiology Investigations. B 29 Issue 7.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887595/B_29i7.pdf
4. Public Health England, 2015. 'Identification of Enterobacteriaceae'. UK Standards for Microbiology Investigations. ID 16 Issue 4.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf

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Bibliografi

1. Foster, Timothy J. 1996. "Staphylococcus." In Medical Microbiology. 4th Edition., 4th ed. Galveston. <https://www.ncbi.nlm.nih.gov/books/NBK8448/>.

2. Public Health England. 2020a. 'Investigation of Staphylococcus species, Micrococcus species and Rothia species.' UK Standards for Microbiology Investigations. ID 7 Issue 4.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887570/UK_SMI_ID_07_i4.pdf
3. Public Health England, 2020b. 'Investigation of specimens for screening for MRSA'. UK Standards for Microbiology Investigations. B 29 Issue 7.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/887595/B_29i7.pdf
4. Public Health England, 2015. 'Identification of Enterobacteriaceae'. UK Standards for Microbiology Investigations. ID 16 Issue 4.
https://assets.publishing.service.gov.uk/government/uploads/system/uploads/attachment_data/file/423601/ID_16i4.pdf

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© 2022 Thermo Fisher Scientific Inc. Med ensamrätt.
Alla övriga varumärken tillhör Thermo Fisher Scientific Inc. och dess dotterbolag.