Chromatic™ E.coli O157

Chromogenic medium for the isolation and differentiation of *Escherichia coli* O157 in clinical and food samples.

DESCRIPTION

Chromatic™ E.coli O157 is a chromogenic medium used for the isolation and differentiation of *Escherichia coli* O157 in clinical specimens and food samples.

TYPICAL FORMULA

	(g/l)
Peptone	15.00
Casein Peptone	5.00
Yeast Extract	5.00
Meat Extract	1.00
Sodium Chloride	5.00
Sodium Deoxycholate	0.50
Chromogenic Mix	2.25
Agar	15.00

Final pH 7.0 ± 0.2 at 25°C

METHOD PRINCIPLE

Peptones, yeast and meat extracts are sources of amino acids, nitrogen, minerals and vitamins. Sodium chloride maintains the osmotic balance of the medium. Sodium deoxycholate inhibits growth of Gram-positive bacteria. The chromogenic mix allows the identification of the microorganisms on the basis of the color and morphology of the colonies. Agar is the solidifying agent.

PREPARATION

Dehydrated medium Suspend 49 g of the powder in 1 liter of distilled or deionized water. Mix well. Heat to boil shaking frequently until completely dissolved. DO NOT AUTOCLAVE. Cool to 45-50°C. In aseptic conditions dispense in Petri dishes. Allow the medium to solidify with lids partially removed.

Medium in bottles Melt the content of the bottle in a water bath at 100°C (loosing the cap partially removed) until completely dissolved. Then screw the cap and check the homogeneity of the dissolved medium, if it is the case turning the bottle upside down. Cool at 45-50°C, mix well avoiding foam formation and aseptically distribute into Petri dishes.

TEST PROCEDURE

Inoculate the plates by streaking directly the sample onto the agar surface. Alternatively, plates may be inoculated from pre-enrichments. Incubate aerobically at 37°C for 18-24 h.

INTERPRETING RESULTS

At the end of incubation, examine for typical mauve or mauve-purple colonies of *Escherichia coli* O157. Gram-positive organisms are completely inhibited. Gram-negative organisms, other than *E. coli* O157, are either inhibited or produce colorless, blue, green, or blue-green colonies. Some *Salmonella* spp. could give false positive results. Final identification should be confirmed by immunoserological tests.

See pictures in Appendix I.

APPEARANCE

Dehydrated medium: free-flowing homogeneous, beige to faint pink.

Bottle: opaque, beige medium.

Prepared medium: slightly opalescent, light beige.

STORAGE

The powder is very hygroscopic, store the powder at 10-30°C, in a dry environment, in its original container tightly closed. Store bottles and prepared plates at 2-8°C away from light. Do not use the product beyond its expiry date on the label or if product shows any evidence of contamination or any sign of deterioration.

SHELF LIFE

Dehydrated medium: 2 years.
 Medium in bottle: 1 year.
 Ready-to-use plates: 4 months.

QUALITY CONTROL

Plates are inoculated with the microbial strains indicated in the QC table.
 Inoculum for productivity: 50-100 CFU.
 Inoculum for selectivity: 10⁴-10⁶ CFU.
 Inoculum for specificity: 10³-10⁴ CFU
 Incubation conditions: aerobically at 35±2°C for 18-24 h.

QC Table.

Microorganism		Growth	Specification
<i>Escherichia coli</i> O157	ATCC® 35150	Good	Mauve or mauve-purple colonies
<i>Escherichia coli</i>	ATCC® 25922	Good	Blue colonies
<i>Klebsiella pneumoniae</i>	ATCC® 13883	Good	Blue colonies
<i>Staphylococcus aureus</i>	ATCC® 25923	Inhibited	---
<i>Enterococcus faecalis</i>	ATCC® 19433	Inhibited	---

WARNING AND PRECAUTIONS

The product does not contain hazardous substances in concentrations exceeding the limits set by current legislation and therefore is not classified as dangerous. It is nevertheless recommended to consult the safety data sheet for its correct use. The product is intended for *In vitro* diagnostic use and must be used only by properly trained operators.

DISPOSAL OF WASTE


Disposal of waste must be carried out according to national and local regulations in force.

BIBLIOGRAPHY

1. L.M. Durso and J.E. Keen USDA, ARS, US Meat Animal Research Center, Clay Center, NE, USA Publication: Journal of Applied Microbiology ISSN 1364-5072: 2007.

PRESENTATION		Contents	Ref.
Chromatic™ E.coli O157	90 mm ready-to-use plates	20 plates	11610
Chromatic™ E.coli O157	Bottles	6 x 200 ml bottles	482190
Chromatic™ E.coli O157	Dehydrated medium	500 g of powder	610614
Chromatic™ E.coli O157	Dehydrated medium	100 g of powder	620614

TABLE OF SYMBOLS

LOT Batch code	IVD <i>In vitro</i> Diagnostic Medical Device	 Manufacturer	 Use by	 Fragile, handle with care
REF Catalogue number	 Temperature limitation	 Contains sufficient for <n> tests	 Caution, consult Instruction For Use	 Do not reuse

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