

---

# Shaker Incubator SKI 4

Item Code: 41102012



## DATA SHEET

### Features

---

Accuracy (stirrer speed)	1 rpm
Power supply	230 V
Capacity	4 L
Heating accuracy control	0,1 °C
Plate dimension	320 x 320 mm
Dimension	500 x 610 x 470 mm
Stir excursion	20 mm
Frequency	50 / 60 Hz Hz
Warranty	2 anni
Image 2	<a href="https://app.onpage.it/api/storage/6baee01c6b26a7cd4a6ec3eccc843d0b7c86">https://app.onpage.it/api/storage/6baee01c6b26a7cd4a6ec3eccc843d0b7c86</a>
Temperature indication	Yes
Time indication	Yes
Speed indication	Yes
Plate material	Metal
Stir movement	Orbital
Name	SKI 4
Regulations	2006/95/EC -2004/108/EC - EN61326-1:2006 - EN61010-1:2010 - EN61010-2-010:2003
Weight	40 kg
Absorbed power	500 W
Heating power	500 W
Over Temperature Protection	Yes
Heating Temperature Range	+ 5 Ta ... 60 °C
Timer	99:59 o continuo hh:mm
Stirring speed	40 ... 300 rpm

---