

---

# HB 120-C Dry bath for laboratory

Item Code: 22010013



## DATA SHEET

### Features

---

Accuracy (Temperature)	±0.5 °C
Power supply	110/220, 50-60Hz V
Alarms	Yes, acoustic/visual
Dimension	290x175x85 mm
Display	LED
Warranty	2 anni
IP rating	21 IP
Image 2	<a href="https://app.onpage.it/api/storage/2526b9c49071a25450dbb64def64ff8d317f10">https://app.onpage.it/api/storage/2526b9c49071a25450dbb64def64ff8d317f10</a>
Name	HB 120
Weight	3 kg
Heating power	160 W
Over Temperature Protection	Yes, 140°C
Calibration points	1
temperature measurement range	+25...120 °C
Working temperature	10...40 °C
Timer	1min-99h:59min hh:mm
Working humidity	< 80% %RH

---