

Technical Datasheet

Article No. 2935

Nitrate – Test strips 0 - 500 mg/l, 100 pieces, semi-quantitative

Single-use test strips for the rapid, semi-quantitative detection of nitrate ions in solutions.

Parameter	Value
Measuring range	0 - 500 mg/l
Gradation	0 - 10 - 25 - 50 - 100 – 250 - 500 mg/l
Colour reaction	White-pale pink to violet
Shelf life	2 years from date of manufacture Use opened packs within 3 months.

Instructions for use

1. Immerse the test strip in the solution for 2 seconds.
2. Shake off excess solution.
3. Wait 30 seconds.
4. Compare the colour produced on the strip with the colour chart.

Storage instructions

Unused strips must be stored in the original packaging. This must be kept tightly closed in a cool, dry place where it is not exposed to chemical vapours.

Version No. 02/2024-11-18