

Z600008

 Partager

Hellma® absorption cuvettes, standard cells, Macro



Special Optical Glass, spectral range 360-2,500 nm, pathlength 10 mm, chamber volume 3,500 μ L

Synonyme(s) : Hellma Macro Absorption Cuvettes

PROPRIÉTÉS

Matériaux	optical glass
Description	2 polished windows
Fabricant/nom de marque	Hellma 6030-OG
Paramètres	360-2,500 nm spectral range
Technique(s)	UV/Vis spectroscopy: suitable photometry: suitable
H x l x P	45 mm x 12.5 mm x 12.5 mm
Volume de la chambre	3,500 μ L
Longueur du trajet optique	10 mm
Méthode de détection	photometric

DESCRIPTION

Description générale

without lid

Conditionnement

Sold as packages containing a single cell (1EA = 1 cell)

Informations légales

Hellma is a registered trademark of Hellma GmbH