

XS Sensor Polymer S7

Item Code: 32200223



DATA SHEET

Features

Connection	Screw head S7
Description	Combined pH electrode, screw head S7. General purpose glass body. Cylindrical membrane. Polymer filling without porous septum. For fast response with high accuracy in samples that may contaminate normal diaphragms. Maintenance-free. Dimensions (WxD) mm 120x12 - pH 0...14 - Temperature - 10...80 °C
Diaphragm	Open hole
Body diameter	12 mm
Electrolyte	Polymer
Image 2	https://app.onpage.it/api/storage/29460e289afe69a65b49e3843d74c8ff0d2aC
Body length	120 mm
Body material	Glass
Name	Polymer
pH measuring range	0 ... 14 pH
temperature measurement range	-10 ... +80 °C
Reference system	Ag / AgCl Encapsulated