

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

## **Product Specification**

Product Name:

Sodium hydroxide - pellets, semiconductor grade, 99.99% trace metals basis

 Product Number:
 306576

 CAS Number:
 1310-73-2

 MDL:
 MFCD00003548

Formula: HNaO

Formula Weight: 40.00 g/mol

NaOH

TEST Specification

Appearance (Color) White

Appearance (Form) Pellets

X-Ray Diffraction Conforms to Structure

Titration by HCL  $\geq$  95.0 %

(% NaOH)

Impurity < 1.0 % % Na2CO3

Trace Metal Analysis ≤ 150.0 ppm

Purity Meets Requirements

99.99% Based on Trace Metals Analysis

Specification: PRD.2.ZQ5.10000004617

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.