



Pyrex™ Screwcap Spare Inserts

For use with Black SVL™ PBT screwcaps

Brand: Pyrex 4513/42M

Product

Réf.: 274323

Description

- PTFE faced sealing disc for chemical resistance
- Autoclavable up to 121°C

Specifications

| | |
|--------------|--------------|
| Closure Size | 42 mm |
| Description | Spare insert |
| Quantity | 10 Pack |

| | |
|--------------|--------------|
| Autoclavable | Autoclavable |
| Color | White |



Distribué par **Humeau**

Z.A. Hauts de Couëron - 11 rue des Imprimeurs
CS 24220 - 44220 Couëron - France
t. : +33 (0)2 40 93 53 53
info@humeau.com - www.humeau.com