

Mini CENTRIFUGE D2012 with Rotor A12-2

Item Code: *22000012*

DATA SHEET

Features

Power supply	220, 50/60 Hz V
Capacity max. (fixed angle rotor)	12 x 2 ml
Consumption	120 W
Dimension	255 x 245 x 140 mm
Standard - directives	2014/35/EU ; 2014/30/EU ; 2015/863/EU ; EN61010-1 ; EN61326-1 ; EN61000-3-2 ; EN 61000-3-3 ; EN 61010-2-020
RCF display	Yes
Spin function	Yes
Warranty	2 anni
Image 2	https://app.onpage.it/api/storage/76098aae14e014a1137a07edf754cf9bd1baca
RCF setting	Yes
RPM settings	Yes
Name	D 2012
Weight	46 kg
RCF max (fixed angle rotor)	15100 x g
RCF max	15100 xg
Noise	56 dB
Control system	Microprocessor and LCD display
Timer	30 sec. - 99 min hh:mm
Type of motor	Brushless induction
maximum speed	15000 rpm