

Product Specifications

CL00.0359**tri-Calcium phosphate v.p.**

tri-Calcio fosfato t.p.
tri-Calcio fosfato e.p.
tri-Calciumfosfaat z.z.
tri-Calciumphosphat Rst.
tri-Calcium phosphate t.p.

Mol weight: 502.32 g/mol
Density: 3,14 g/cm³
HS Nr: 28352600
CAS Nr: 7758-87-4
UN:
ADR:

For laboratory use, BP, E 341
Ca₅(PO₄)₃OH - Ca₃(PO₄)₂ > 90%

Specifications:

Calcium Phosphate	>90%
Insoluble Matter (Non Solubles)	<0.2%
Water	<2.5%
Loss on ignition	<2.5%
Arsenic (As)	<0.0002%
Copper (Cu)	<0.002%
Iron (Fe)	<0.04%
Lead (Pb)	<0.0005%
Heavy Metals as Lead (Pb)	<0.003%
Chloride (Cl)	<0.3%
Fluoride (F)	<0.005%
Sulfate (SO ₄)	<0.5%