

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Research[®] plus Pipettes

Eppendorf Research [®] plus/ Eppendorf Research [®] plus Move It [®]	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	0.1 – 20 µL, light pink	5 – 100 µL, light yellow	2 – 20 µL, yellow	2 – 100 µL, yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.25 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise	GEloader [®] / Microloader™ 0.5-20 µL
epT.I.P.S. [®]	■		■	■	■				■		■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] LoRetention [®]	■			■	■				■		■	■	■	■	■	■	■	■	■
ep Dualfilter T.I.P.S. [®]	■	■		■	■	■	■	■	■	■	■	■	■	■	■	■	■	■	■
ep Dualfilter T.I.P.S. [®] LoRetention [®]	■			■	■				■		■	■	■	■	■	■	■	■	■
ep Dualfilter T.I.P.S. [®] SealMax [®]	■			■	■				■		■	■	■	■	■	■	■	■	■
Fixed volume, single-channel pipettes																			
10 µL	•	■	■	•															•
20 µL			•	■															•
10 µL	■																		
20 µL	■						■	■	■	•									
25 µL/50 µL/100 µL	■						■	■	■	•									
200 µL	■							■	■	•									
200 µL/250 µL/500 µL/1,000 µL	■								■	•	•	•							
Variable volume, single-, 8- and 12-channel pipettes																			
0.1 – 2.5 µL	■	•	•	•															
0.5 – 10 µL	■	•	■	•															•
2 – 20 µL	■	•	•	■															•
2 – 20 µL	■						■	■	■	•									
10 – 100 µL	■						•	■	■	•									
20 – 200 µL	■						•	■	■	•									
30 – 300 µL	■						•	■	■	■									
100 – 1,000 µL	■										■	•	•						
120 – 1,200 µL	■										•	•	■						
0.25 – 2.5 mL	■													■					
0.5 – 5 mL															■				
1 – 10 mL	■															•	■	•	
Variable volume, 16- and 24-channel pipettes																			
1 – 20 µL	■				■														
5 – 100 µL	■					■													
Move It[®], variable volume, 4-, 6-, 8- and 12-channel pipettes																			
1 – 20 µL	■				■														
5 – 100 µL	■					■													
30 – 300 µL	■						•	•	•	■									
120 – 1,200 µL	■										•	•	■						

■ Test tip (other tips may require adjustment) • Full volume range ◊ Limited volume range



Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

For more information go to
www.eppendorf.com/eptips

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Reference[®] 2 Pipettes

Eppendorf Reference [®] 2	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	2 – 20 µL, yellow	2 – 100 µL, yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.25 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise	GEloader [®] / Microloader™ 0.5 – 20 µL
epT.I.P.S. [®]	■		■	■			■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] LoRetention [®]	■			■			■		■		■		■		■		■
ep Dualfilter T.I.P.S. [®]	■	■		■	■	■	■	■	■		■	■	■	■	■		
ep Dualfilter T.I.P.S. [®] LoRetention [®]	■			■		■		■	■								
ep Dualfilter T.I.P.S. [®] SealMax [®]	■			■		■		■	■								
Fixed volume, single-channel pipettes																	
1 µL/2 µL	■	•	•														
5 µL/10 µL	■	•	■	•													•
20 µL	■		•	■													•
10 µL	■				■	■	■	•									
20 µL	■				■	■	■	•									
25 µL/50 µL/100 µL	■					■	■	•									
200 µL	■						■	•									
200 µL/250 µL/500 µL/1,000 µL	■								■	•	•						
2 mL/2.5 mL	■											■					
Variable volume, single-, 8- and 12-channel pipettes																	
0.1 – 2.5 µL	■	•	•	•													•
0.5 – 10 µL	■	•	■	•													•
2 – 20 µL	■	•	•	■													•
2 – 20 µL	■				◆	◆	■	•									
10 – 100 µL	■				◆	■	■	•									
20 – 200 µL	■				◆	■	■	•									
30 – 300 µL	■				◆	◆	◆	■									
100 – 1,000 µL	■								■	•	•						
0.25 – 2.5 mL	■											■					
0.5 – 5 mL	■												■	•			
1 – 10 mL	■													■	•	■	•

■ Test tip (other tips may require adjustment) • Full volume range ◆ Limited volume range



Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

For more information go to
www.eppendorf.com/eptips

Combination of epT.I.P.S.[®] Pipette Tips with Eppendorf Xplorer[®]/Xplorer[®] plus Pipettes

Eppendorf Xplorer [®] /Eppendorf Xplorer [®] plus/ Eppendorf Xplorer [®] plus Move It [®]	0.1 – 10 µL, dark gray	0.1 – 10 µL M, medium gray	0.1 – 20 µL, medium gray	0.5 – 20 µL L, light gray	0.1 – 20 µL, light pink	5 – 100 µL, light yellow	2 – 20 µL, yellow	2 – 100 µL, yellow	2 – 200 µL, yellow	20 – 300 µL, orange	50 – 1,000 µL, blue	50 – 1,250 µL, green	50 – 1,250 µL L, dark green	0.25 – 2.5 mL, red	0.1 – 5 mL, violet	0.2 – 5 mL L, violet	0.5 – 10 mL, turquoise	0.5 – 10 mL L, turquoise	GEloader [®] / Microloader [™] 0.5-20 µL
epT.I.P.S. [®]	■		■	■	■				■	■	■	■	■	■	■	■	■	■	■
epT.I.P.S. [®] LoRetention [®]	■			■					■		■			■		■			■
ep Dualfilter T.I.P.S. [®]	■	■		■	■	■	■	■	■	■	■		■	■	■	■	■	■	
ep Dualfilter T.I.P.S. [®] LoRetention [®]	■			■				■	■	■	■								
ep Dualfilter T.I.P.S. [®] SealMax [®]	■			■				■	■	■	■								
Variable volume, single-, 8- and 12-channel pipettes																			
0.5 – 10 µL	■	●	■	●															●
1 – 20 µL	◆	◆	●	■															●
5 – 100 µL							◆	■	■	●									
10 – 200 µL							◆	◆	■	●									
15 – 300 µL							◆	◆	◆	■									
50 – 1,000 µL											■	●	●						
50 – 1,200 µL											◆	■	●						
0.1 – 2.5 mL														■					
0.2 – 5 mL															■	●			
0.5 – 10 mL																	■	●	
Variable volume, 16- and 24-channel pipettes																			
1 – 20 µL					■														
5 – 100 µL						■													
Move It[®], variable volume, 4-, 6-, 8- and 12-channel pipettes																			
1 – 20 µL					■														
5 – 100 µL						■													
15 – 300 µL							◆	◆	◆	■									
50 – 1,200 µL											◆	■	●						

■ Test tip (other tips may require adjustment) ● Full volume range ◆ Limited volume range



Your local distributor: www.eppendorf.com/contact
 Eppendorf SE · 22339 Hamburg · Germany
eppendorf@eppendorf.com · www.eppendorf.com

www.eppendorf.com

For more information go to
www.eppendorf.com/eptips