

# PX Vio Simulator - KIT

Item Code: 50111312



Accuracy (pH)	±0,02 pH
Power supply	4 batteries AA 1,5 V V
Battery supply	Yes
Measuring range	pH: -2,00 ... 16,00 (0,01 pH) mV: -1999 ... +1999 (1 mV)
Dimension	185 x 85 x 45 mm
Display	LCD high definition, backlight
Supplied as standard	1 meter BNC/BNC cable, coloured pH buffer solutions, operating instructions, test report and carrying case.
IP rating	57 IP
Image 2	<a href="https://app.onpage.it/api/storage/8f722f2fa5fd91901568320aac9c2d6d06bd6c">https://app.onpage.it/api/storage/8f722f2fa5fd91901568320aac9c2d6d06bd6c</a>
Inputs	1 for BNC
Regulations	CE
Weight	0,4 kg
Measuring range mV	-1999 ... +1999 mV
pH measuring range	-2,00 ... 16,00 pH
pH resolution	0,01 pH
Auto off	20 minutes (can be deactivated)
Working temperature	0 ... +60 °C