
1X500ml XS Solution Cond. 12880

$\mu\text{S}/\text{cm} \pm 1\%/25^\circ\text{C}$, with certificate traceable to NIST

Item Code: 51100643



DATA SHEET

Features

Accuracy	$\pm 1\%$
Certificate	Yes
Colour	Colourless
Composition	Potassium chloride in aqueous solution
Pack	500 ml with integrated calibration compartment
Description	11670 $\mu\text{S}/\text{cm}$ (20°C) - N.I.S.T. Accuracy $\pm 1\%$
Name	EC 12,880 $\mu\text{S}/\text{cm}$
Nominal value	12880 $\mu\text{S}/\text{cm}$ (25 °C)
Volume	500 ml
