
1X500ml XS Professional 1413 μ S/

cm

$\pm 1\%$ / 25°C, with certificate

DFM/DANAK

Item Code: 51300343



DATA SHEET

Features

Accuracy	$\pm 1\%$
Certificate	Yes
Colour	Colourless
Composition	Potassium chloride in aqueous solution
Pack	500 ml with integrated calibration compartment
Description	1278 μ S/cm (20°C) - DFM. Accuracy $\pm 1\%$
Name	EC 1413 μ S/cm
Nominal value	1413 μ S/cm (25 °C)
Volume	500 ml
