

Safety Data Sheet

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product name : KWIK-STIK™ Swab Hydrating Fluid

Trade names : KWIK-STIK™ KWIK-STIK™ F

KWIK-STIK™ Plus Lab-Elite™ CRM

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Use of the substance/mixture : Hydrating fluid

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Microbiologics, Inc.

200 Cooper Avenue North

Saint Cloud, MN 56303

+1.320.253.1640

1.4. Emergency telephone number

24 hour Emergency Number: +44 1865 407333 (Carechem)

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substance

Not applicable

3.2. Mixture

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Water	(CAS No) 7732-18-5 (EC no) 231-791-2	> 99	Not classified
Sodium chloride	(CAS No) 7647-14-5 (EC no) 231-598-3	< 1	Not classified
Sodium phosphate dibasic	(CAS No) 7558-79-4 (EC no) 231-448-7	< 1	Not classified
Sodium thioglycolate	(CAS No) 367-51-1 (EC no) 206-696-4	< 1	Acute Tox. 4 (Oral), H302 Skin Irrit. 2, H315 Eye Irrit. 2, H319
Phosphoric acid, potassium salt (1:1)	(CAS No) 7778-77-0 (EC no) 231-913-4	< 1	Not classified
Magnesium chloride	(CAS No) 7786-30-3 (EC no) 232-094-6	< 1	Not classified
Calcium chloride	(CAS No) 10035-04-8 (EC no) 600-075-5	< 1	Acute Tox. 4 (Oral), H302 Eye Irrit. 2, H319 STOT SE 3, H335

Full text of R- and H-statements: see section 16

Safety Data Sheet

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : No specific first aid measures needed for this route of exposure.

First-aid measures after skin contact : Wash with soap and water. Seek medical assistance if irritation develops or persists.

First-aid measures after eye contact : Rinse cautiously with water for several minutes. Remove contact lenses, if present and easy to

do. Continue rinsing. If irritation persists, get medical advice/attention.

First-aid measures after ingestion : Avoid hand to mouth contact. If ingested, seek medical advice.

4.2. Most important symptoms and effects (acute and delayed)

Symptoms/injuries after inhalation : None anticipated under normal product handling conditions.

Symptoms/injuries after skin contact : May cause irritation. Symptoms/injuries after eye contact : May cause irritation.

Symptoms/injuries after ingestion : May be harmful if swallowed.

4.3. Immediate medical attention and special treatment, if necessary

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Use suitable extinguishing media for surrounding fire.

Unsuitable extinguishing media : None.

5.2. Special hazards arising from the substance or mixture

Fire hazard : None known. Explosion hazard : None known.

5.3. Advice for firefighters

Protection during firefighting : Firefighters should wear full protective gear.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

No additional information available

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

For containment : Stop the flow of material, if this is without risk.

Methods for cleaning up : If hydration of the lyophilized microorganism preparation has not occured, no action is required.

If hydration has occured, please see LIT.115 Biohazard Cleanup on our website at

www.microbiologics.com.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling

: The hydrating fluid is a sterile fluid and, by itself, does not pose any hazardous threats. When used to hydrate the lyophilized microorganism preparation, the hydrating fluid will create a suspension that does contain microorganisms, which under certain conditions, could lead to an infectious process.

Proper techniques must be employed to avoid exposure and contact with microorganism growth, and rehydrated pellet suspensions. The microbiology laboratory must be equipped, and have the facilities to receive, process, maintain, store and dispose of biohazard material. The microbiology laboratory personnel using these devices must be trained, experienced, and demonstrate proficiency in processing, maintaining, storing and disposing of biohazard material.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store the product at 2°C - 8°C in the original sealed container.

7.3. Specific end use(s)

No additional information available

SDS.2150.ENG.EU Rev B English Page 2/6 2019.JULY.11

Safety Data Sheet

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Sodium chloride (7647-14-5)		
Latvia	OEL TWA (mg/m³)	5 mg/m³
Lithuania	IPRV (mg/m³)	5 mg/m³

8.2. Exposure controls

Appropriate engineering controls : Local exhaust and general ventilation must be adequate to meet exposure standards. Good

laboratory practices must be observed and followed.

Hand protection : Wear general protective gloves.

Eye protection : Safety glasses with side shields.

Skin and body protection : Wear moisture impervious aprons and safety footwear.

Respiratory protection : When undertaking procedures that are likely to give rise to infectious aerosols, a Class 1

microbiological biological safety cabinet should be used.

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Liquid. Each unit contains a reservoir of hydrating fluid in the stick.

Colour : Colourless
Odour : Odourless

No data available Odour threshold No data available pΗ Relative evaporation rate (butylacetate=1) : No data available Melting point : No data available No data available Freezing point Boiling point No data available : No data available Flash point Auto-ignition temperature : No data available Decomposition temperature No data available Flammability (solid, gas) No data available Vapour pressure : No data available : No data available : No data available No data available : No data available

Relative vapour density at 20 °C : No data available Relative density : No data available Solubility : No data available Log Pow : No data available Viscosity, kinematic : No data available Viscosity, dynamic : No data available Explosive properties : No data available Oxidising properties : No data available Explosive limits : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

Stable under normal ambient and anticipated storage and handling conditions.

10.3. Possibility of hazardous reactions

Will not occur.

10.4. Conditions to avoid

None.

10.5. Incompatible materials

Not determined.

10.6. Hazardous decomposition products

Not determined.

SDS.2150.ENG.EU Rev B English Page 3/6 2019.JULY.11

Safety Data Sheet

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity (oral) : Not classified
Acute toxicity (dermal) : Not classified
Acute toxicity (inhalation) : Not classified

Water (7732-18-5)

LD50 oral rat > 90 ml/kg

Sodium chloride (7647-14-5)

LD50 oral rat	3 g/kg

LC50 inhalation rat (mg/l) > 42 g/m³ (Exposure time: 1 h)

Magnesium chloride (7786-30-3)

LD50 oral rat 2800 mg/kg

Sodium phosphate dibasic (7558-79-4)

LD50 oral rat 17 g/kg

Phosphoric acid, potassium salt (1:1) (7778-77-0)

LD50 oral rat 3200 mg/kg
Skin corrosion/irritation : Not classified

Serious eye damage/irritation : Not classified
Respiratory or skin sensitisation : Not classified
Germ cell mutagenicity : Not classified
Carcinogenicity : Not classified
Reproductive toxicity : Not classified

Specific target organ toxicity (single exposure) : Not classified

Specific target organ toxicity (repeated

exposure)

: Not classified

Aspiration hazard : Not classified

SECTION 12: Ecological information

12.1. Toxicity

Acute aquatic toxicity : Not classified Chronic aquatic toxicity : Not classified

Sodium chloride (7647-14-5)	
LC50 fish 1	5560 - 6080 mg/l (Exposure time: 96 h - Species: Lepomis macrochirus [flow-through])
LC50 fish 2	12946 mg/l (Exposure time: 96 h - Species: Lepomis macrochirus [static])
EC50 Daphnia 1	1000 mg/l (Exposure time: 48 h - Species: Daphnia magna)
EC50 Daphnia 2	340.7 - 469.2 mg/l (Exposure time: 48 h - Species: Daphnia magna [Static])

	LOJU Daprillia Z	540.7 - 409.2 mg/1 (Exposure time: 40 ft - Species: Daprillia magna [Static])
Magnesium chloride (7786-30-3)		
	LC50 fish 1	1970 - 3880 mg/l (Exposure time: 96 h - Species: Pimephales promelas [static])
	EC50 Daphnia 1	140 mg/l (Exposure time: 48 h - Species: Daphnia magna [Static])
	EC50 72h algae (1)	2200 mg/l (Species: Desmodesmus subspicatus)

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

Sodium chloride (7647-14-5)	
BCF fish 1	(no bioaccumulation)

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

No additional information available

SDS.2150.ENG.EU Rev B English Page 4/6 2019.JULY.11

Safety Data Sheet

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Product/Packaging disposal recommendations : Dispose of contents/container in accordance with local/regional/national/international

regulations.

SECTION 14: Transport information

In accordance with ADR / RID / IMDG / IATA / ADN

14.1. UN number

UN-No. (ADR) : Not applicable
UN-No. (IMDG) : Not applicable
UN-No. (IATA) : Not applicable
UN-No. (ADN) : Not applicable
UN-No. (RID) : Not applicable

14.2. UN proper shipping name

Proper Shipping Name (ADR) : Not applicable
Proper Shipping Name (IMDG) : Not applicable
Proper Shipping Name (IATA) : Not applicable
Proper Shipping Name (ADN) : Not applicable
Proper Shipping Name (RID) : Not applicable

14.3. Transport hazard class(es)

ADR

Transport hazard class(es) (ADR) : Not applicable

IMDG

Transport hazard class(es) (IMDG) : Not applicable

IATA

Transport hazard class(es) (IATA) : Not applicable

ADN

Transport hazard class(es) (ADN) : Not applicable

RID

Transport hazard class(es) (RID) : Not applicable

14.4. Packing group

Packing group (ADR) : Not applicable
Packing group (IMDG) : Not applicable
Packing group (IATA) : Not applicable
Packing group (ADN) : Not applicable
Packing group (RID) : Not applicable

14.5. Environmental hazards

Dangerous for the environment : No Marine pollutant : No

Other information : No supplementary information available

14.6. Special precautions for user

- Overland transport

Not applicable

- Transport by sea

Not applicable

- Air transport

Not applicable

- Inland waterway transport

Not applicable

SDS.2150.ENG.EU Rev B English Page 5/6 2019.JULY.11

Safety Data Sheet

- Rail transport

Not applicable

14.7. Transport in bulk according to Annex II of MARPOL and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no substances with Annex XVII restrictions Contains no substance on the REACH candidate list

Contains no REACH Annex XIV substances

15.1.2. National regulations

Germany

AwSV/VwVwS Annex reference : Water hazard class (WGK) 3, severe hazard to waters (Classification according to AwSV,

Annex 1)

12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende stoffen : Magne

SZW-lijst van mutagene stoffen

NIET-limitatieve lijst van voor de voortplanting

giftige stoffen – Borstvoeding

NIET-limitatieve lijst van voor de voortplanting

giftige stoffen - Vruchtbaarheid

NIET-limitatieve lijst van voor de voortplanting

giftige stoffen - Ontwikkeling

: Magnesium chloride is listed: Magnesium chloride is listed

Magricolam ornorae io notea

: None of the components are listed

: None of the components are listed

: None of the components are listed

15.2. Chemical safety assessment

No additional information available

SECTION 16: Other information

Full text of R-, H- and EUH-statements:

Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Skin Irrit. 2	Skin corrosion/irritation, Category 2
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H302	Harmful if swallowed
H315	Causes skin irritation
H319	Causes serious eye irritation
H335	May cause respiratory irritation

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product