



Product Specifications

CL02.2609

Sulfuric acid 90-91 weight % solution

Schwefelsäure 90-91 w/w % Lösung
Acide sulfurique 90-91 poids % solution
Acido solforico 90-91 p/p % soluzione
Zwavelzuur 90-91 gew. % oplossing
Acido sulfúrico 90-91 p/p % solución

**For laboratory use / milk analysis according GERBER
1633 g H₂SO₄ / l H₂O**

Mol weight: 98.08 g/mol
Density: 1,81 g/ml
HS Nr: 28070000
CAS Nr: 7664-93-9
UN: 1830
ADR: 8,II



Specifications:

Colour < 10 APHA