



# Product Specifications

**CL05.2622**

TITREX

**Sulfuric acid 5 mol/l**

Acido solforico 5 mol/l  
Zwavelzuur 5 mol/l  
Acido sulfúrico 5 mol/l  
Acide sulfurique 5 mol/l  
Schwefelsäure 5 mol/l

**For laboratory use, Titrations**  
**490.4 g H<sub>2</sub>SO<sub>4</sub> / l H<sub>2</sub>O = 10 N (±0.02/20°C)**  
**Traceability: Reference standards acc. NIST SRM 84I**

Mol weight: 98.08 g/mol  
Density: 1,29 g/ml  
HS Nr: 28070000  
CAS Nr: 7664-93-9  
UN: 2796  
ADR: 8,II



## **Specifications:**