

Product Details

This product is supplied with a **certificate of analysis**, and is produced in a laboratory with the following accreditations:
ISO 9001, ISO 17025 & ISO Guide 34

Category: **EUROPEAN PHARMACOPOEIA REAGENT**

Family: **EUROPEAN PHARMACOPOEIA**

Ref: **RP 035 100**

SHIP WEIGHT:0,125 Kg

Description: **Bromine Water.**

Eur Ph Ref.: **1012402**

Volume (ml): **100ml Bottle(s)**

Shelf life: **3 Months.**

Specification from Chapter 4.1.1. of the European Pharmacopoeia: **Bromine water. 1012402. Shake 3 ml of bromine R with 100 ml of water R to saturation. Storage: over an excess of bromine R, protected from light.**

SDS Name: **Bromine Water. European Pharmacopoeia Ref: 1012402 (RP 035)**

Pictograms: **GHS05 ; GHS07**

Signal Word: **Danger**

H Phrases: **H315 ; H318 ; H332**

P Phrases: **P261;P264;P271;P280;P302+P352;P304+P340;P305+P351+P338;P310;P332
+P313;P362+P364**

EUH Statements: **N/A**

UN Number: **N/A**

Packing Group: **N/A**

Chemical Class ADR: **N/A**

ADR Labels: **N/A**

Proper Shipping Name: **N/A**

PSN Addition: **N/A**

DGN Description: **N/A**