



## ZPACK13415030 GLASS BOTTLE LIQUID 250 ML DIN 32



### SPECIFICATIONS

Material :	Glass
Neck Type :	DIN 32
Filling Opening :	32,00 mm
Height :	155,30 mm
Length :	63,00 mm
Wide :	63,00 mm

### AVAILABLE ACCESSORIES

Adaptator for dispensor S40 to GL 45 :	217-9271
Key for opening :	84544.001

Remark : Packaging materials and specifications shown may be modified without preceding notice  
Version dd. 22 mai 2012

tel/tél +32(0)16 385 011  
fax +32(0)16 385 385  
[info@be.vwr.com](mailto:info@be.vwr.com)  
[www.vwr.com](http://www.vwr.com)

Researchpark Haasrode 2020  
Geldenaaksebaan 464 B-3001 Leuven  
RPR/RPM Leuven  
BTW/TVA BE 0403.593.343