

3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.com
Email USA: techserv@sial.com
Outside USA: eurtechserv@sial.com

Product Specification

Product Name:

Formula Weight:

BES - ≥99.0% (titration)

Product Number: B9879

CAS Number: 10191-18-1 MFCD00007533 Formula: C6H15NO5S

213.25 g/mol

HO SO₃H

TEST

Appearance (Color)

Appearance (Form)
Solubility (Color)
Solubility (Turbidity)

25 g plus 50 ml of H2O

Infrared spectrum
Water (by Karl Fischer)
A290 UV absorbance

33% W/W

Heavy Metals (as Lead) Recommended Retest Period

5 years

Purity (Titr. by NAOH) anhyd.

Specification: PRD.0.ZQ5.10000023619

Specification

White

Crystalline Powder

Colorless Clear

Conforms to Structure

≤ 0.5 %< 0.05

< = 5 ppm

> 99.0 %

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.