



Sensor type	temperature probe
Order number	51300164
Specifications	
pH range	
zero point (pH)	-
temperature range	-30...100 °C
temperature probe	NTC 30 KΩ
shaft material	stainless steel
membrane glass	-
membrane resistance	-
reference system	-
reference diaphragm	-
reference electrolyte	-
bridge electrolyte	-
electrode head	-
cable length	1.0 m
connector	Cinch
cell constant (cm ⁻¹)	-
cell type	-
indicating metal	-
Storage	
Dry	