


3050 Spruce Street, Saint Louis, MO 63103, USA

Website: www.sigmaaldrich.comEmail USA: techserv@sial.comOutside USA: eurtechserv@sial.com

Product Specification

Product Name:
Titanium(IV) isopropoxide - 97%

Product Number: 205273
CAS Number: 546-68-9
MDL: MFCD00008871
Formula: C₁₂H₂₈O₄Ti
Formula Weight: 284.22 g/mol



TEST

Specification

Infrared spectrum	Conforms to Structure
Specific Gravity @ 25/25 C	0.950 - 0.960
Titanium as TiO ₂	27.5 - 28.3 %
Chloride (Cl)	≤ 65 ppm
Freezing point At Time of Manufacturing	17.5 - 19.5 °C
Color Test At Time of Manufacturing	≤ 25 APHA
Product of Supplier Dorf Ketal	Confirmed

Specification: PRD.0.ZQ5.10000025531

Sigma-Aldrich warrants, that at the time of the quality release or subsequent retest date this product conformed to the information contained in this publication. The current Specification sheet may be available at Sigma-Aldrich.com. For further inquiries, please contact Technical Service. Purchaser must determine the suitability of the product for its particular use. See reverse side of invoice or packing slip for additional terms and conditions of sale.